

How do I visualize my cells in a Thermo Scientific Nunc Cell Factory System?

Multi-layered Thermo Scientific™ Nunc™ Cell Factory™ systems are a proven solution for the cultivation of adherent cells.

An inverted microscope is the perfect solution to monitor the condition of your culture.

To ensure success, the microscope must have a flat stage. When viewing larger formats such as a 10/13 or 40/52 layer, a removable light source is necessary.


Process:



Begin by removing all stage accessories to provide a flat surface for the Cell Factory system to rest upon. If additional space is required, remove the light source.




Place the Cell Factory system on the stage with the desired location to be visualized over the objective lens. If necessary, replace the light source or place an alternative light source above the Cell Factory system.




If the microscope stage is not capable of supporting the weight of the Cell Factory system, it is advisable to use an adjustable stand to carry some of the weight of the Cell Factory system.

Photographs of cells cultured in a Cell Factory System viewed with an inverted microscope.




95% confluent MRC-5 cells

MRC-5 cells are a human diploid cell line commonly used in the manufacturing of vaccines such as Hepatitis A, rubella, varicella and rabies.



80-90% confluent VERO cells

VERO cells are derived from the kidney of an African green monkey and are commonly used for the production of rotavirus, smallpox and inactivated poliovirus.



**60-70% confluent CEF
(Chicken Embryo Fibroblast)**

CEF's are primary cells used in the production of vaccines such as measles, mumps, rabies and influenza. It is also commonly used in the manufacturing of avian vaccines.

www.thermoscientific.com/cellfactory

© 2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

Australia +61 39757 4300 **Austria** +43 1 801 40 0 **Belgium** +32 53 73 42 41 **China** +800 810 5118 or +400 650 5118
France +33 2 2803 2180 **Germany national toll free** 0800 1 536 376 **Germany international** +49 6184 90 6000
India toll free 1800 22 8374 **India** +91 22 6716 2200 **Italy** +39 02 95059 552 **Japan** +81 3 5826 1616
Netherlands +31 76 579 55 55 **New Zealand** +64 9 980 6700 **Nordic/Baltic/CIS countries** +358 9 3291 0200
Russia +7 812 703 42 15 **Spain/Portugal** +34 93 223 09 18 **Switzerland** +41 44 454 12 22
UK/Ireland +44 870 609 9203 **USA/Canada** +1 866 984 3766 **Other Asian countries** +852 2885 4613
Countries not listed +49 6184 90 6000

FAQLSPCFVISUALIZE 0216

Thermo
SCIENTIFIC

A Thermo Fisher Scientific Brand