

Personal Protective Equipment

Hand Protection Disposable



Gloves, polyisoprene, sterile,
hand specific, EN 455
Class 1 Medical



S3



Made from synthetic polyisoprene and hand specific for maximum comfort. Micro roughened surface for optimum grip. Gamma irradiated, AQL 1.0, meets ASTM F1671-97 Viral Penetration Test.

Catalogue No	Size	Pack qty
SAR-681-020Q	5.5	40
SAR-681-030N	6	40
SAR-681-040K	6.5	40
SAR-681-050H	7	40
SAR-681-060E	7.5	40
SAR-681-070B	8	40
SAR-681-080V	8.5	40
SAR-681-090S	9	40

Gloves, polychloroprene, Sterile, Class 100,
ISO Class 5, hand specific,
PPE Category III



S3



The ideal choice for handling cytostatics.

The DermaShield® offers a massive reduction of potential skin reactions and allergies. It minimises the risks of Type IV allergies because it is free of accelerators, and Type I allergies because it is free of natural rubber. Its polychloroprene composition offers a high level of resistance to ethanol and IPA.



Catalogue No	Size	Pack qty
SAR-335-020F	6	200
SAR-335-030C	6.5	200
SAR-335-040W	7	200
SAR-335-050T	7.5	200
SAR-335-060Q	8	200
SAR-335-070N	8.5	200
SAR-335-080K	9	200

Gloves, nitrile, disposable, purple,
PPE Category III

46



- 100% nitrile
- Latex and protein free
- Non-detectable residual accelerators
- Textured fingertips for improved grip
- Powder free
- 0.10mm thickness
- Exceeds the requirements of EN 374 PPE Category III
- Conforms to EN 455 parts 1 & 2
- Ergonomically designed
- Three times more puncture resistant than latex

Nitrile gloves, 250mm length

Catalogue No	Size	Pack qty
FB69260	XX Small	100
FB69261	X Small	100
FB69262	Small	100
FB69263	Medium	100
FB69264	Large	100
FB69265	X Large	100

Extra protection nitrile gloves, 300mm length

Catalogue No	Size	Pack qty
FB69290	XX Small	50
FB69291	X Small	50
FB69292	Small	50
FB69293	Medium	50
FB69294	Large	50
FB69295	X Large	50