



**Fisherbrand**<sup>®</sup>  
QUALITY. RELIABILITY. VALUE.

Focus on pH

*pH-Fix non-bleeding test strips*

# Focus on pH



**Just dip ... and read!**  
**Rapid and reliable pH testing**

 **Fisher  
Scientific**

A Thermo Fisher Scientific Brand



- 85mm x 6mm sticks with different indicator papers sealed onto each test strip
- Provides outstanding measurements safely over the designated pH range
- Colour scale on the test strip box allows exact comparison with the colour developing on the test strip pad
- Indicator is covalently bonded on the pads and does not bleed into the sample – can be left in sample for an indefinite time without contamination
- Hydrophobic plastic handle protects user from sample and allows pH testing in narrow and deep vessels
- Available in a flat box so they fit into every pocket
- An economic and easy standard in a variety of environments such as:
  - Chemical laboratories
  - Industry
  - Dairies
  - Production
  - Universities
  - Schools
  - Food processing

#### Specifications

Shelf life: four years after production · Storage: 4-30°C

Cat. No	pH range	Pack qty
10642751	0.0 to 14.0	100
10333501	4.5 to 10.0	100
10018080	0.0 to 6.0	100
10271751	7.0 to 14.0	100
10353641	3.6 to 6.1	100
11316454	2.0 to 9.0	100
11396454	7.5 to 9.5	100
11386454	6.0 to 7.7	100
11376454	5.1 to 7.2	100
11306464	7.9 to 9.8	100
11346454	0.3 to 2.3	100
10017950	1.7 to 3.8	100