

Balances, analytical, Extend series



Built in application programs: calculation by a factor (multiplication and division), toggling among up to four different units of measure, dynamic weighing, animal weighing (weigh averaging), % weighing, net total formulation, parts counting, totalising.

- New operating concept: easy configuration of the individual settings guided by short, menu driven text prompts (choice of five languages) and simple navigation using cursor keys
- White, backlit, high contrast display with 15mm digits; exceptionally easy to read under any room lighting conditions
- Large, sturdy and nearly frameless, all glass draft shield chamber on 0.1mg models, for optimal viewing of the weighing chamber
- Draft shield doors and all stainless steel parts of the weighing chamber base are easy to remove in just a few steps for cleaning
- Built in, motorised, calibration weight; maximum weighing accuracy and maximum ease of operation at the touch of a key are standard on all analytical balances and on the precision balances ED and CW
- RS232C bidirectional interface port; optional Sartorius cable for connection of the balance to a USB port
- The most advanced microprocessor technology; more performance power for even faster and more stable readings
- ISO/GLP compliant recording of calibration, adjustment and weight values
- All keys have positive click action on a clearly structured control panel
- Digital filter levels of the newest generation offer optimal prerequisites for adapting the balance to the conditions at the place of use

Technical Specification - General

Calibration	Internal
Electrical supply	230V 50/60Hz

Catalogue No	Legacy No	Alt. No	Capacity, g	Readability, mg	Pan diameter, mm
10297352	-	ED224S	220	0.1	90
10790394	-	ED124S	120	0.1	90

Accessory

Catalogue No	Legacy No	Alt. No	Description
10440253	-	6965542	Disposable sample pans

Service you expect

Our product range is only part of the story. You can be assured that you will receive service second to none...

- Large dynamic warehouse means high product availability
- State-of-the-art inventory technology ensures efficient order processing
- Customer Service team dedicated to making your contact with Fisher Scientific easy and efficient
- Knowledgeable field sales specialists to advise you in product choice
- Promotions and offers ensure you receive the best value for money

Contact your local Customer Service team

