

Baths and circulators, water, thermostatic immersion systems, Optima



Grant



T100 thermostat, temperature range 0°C to 100°C

Catalogue No	Legacy No	Alt. No	Capacity, L	Bath material	Dimensions, internal [w x d x h], mm	Dimensions, external [w x d x h], mm
12727809	-	T100-P5	5	Plastic	160 x 240 x 150	200 x 390 x 380
12737809	-	T100-P12	12	Plastic	280 x 325 x 150	350 x 415 x 380
12747809	-	T100-P18	18	Plastic	290 x 510 x 150	350 x 600 x 380
12757809	-	T100-ST5	5	Stainless steel	150 x 300 x 105	180 x 330 x 395
12767809	-	T100-ST12	12	Stainless steel	300 x 325 x 150	330 x 360 x 395
12777809	-	T100-ST18	18	Stainless steel	300 x 505 x 105	330 x 540 x 395
12787809	-	T100-ST26	26	Stainless steel	300 x 505 x 200	330 x 540 x 395
12797809	-	T100-ST38	38	Stainless steel	300 x 690 x 200	333 x 730 x 450

TC120 thermostat, temperature range -20°C to 120°C

Catalogue No	Legacy No	Alt. No	Capacity, L	Bath material	Dimensions, internal [w x d x h], mm	Dimensions, external [w x d x h], mm
12767819	-	TC120-P5	5	Plastic	160 x 240 x 150	200 x 390 x 380
12777819	-	TC120-P12	12	Plastic	280 x 325 x 150	350 x 415 x 380
12787819	-	TC120-P18	18	Plastic	290 x 510 x 150	350 x 600 x 380
12327637	-	TC120-ST5	5	Stainless steel	150 x 300 x 105	180 x 330 x 395
12707829	-	TC120-ST12	12	Stainless steel	300 x 325 x 150	330 x 360 x 395
12717829	-	TC120-ST18	18	Stainless steel	300 x 505 x 105	330 x 540 x 395
12727829	-	TC120-ST26	26	Stainless steel	300 x 505 x 200	330 x 540 x 395
12737829	-	TC120-ST38	38	Stainless steel	300 x 690 x 200	333 x 730 x 450

TX150 thermostat, temperature range -50°C to 150°C

Catalogue No	Legacy No	Alt. No	Capacity, L	Bath material	Dimensions, internal [w x d x h], mm	Dimensions, external [w x d x h], mm
12707839	-	TX150-P5	5	Plastic	160 x 240 x 150	200 x 390 x 380
12717839	-	TX150-P12	12	Plastic	280 x 325 x 150	350 x 415 x 380
12727839	-	TX150-P18	18	Plastic	290 x 510 x 150	350 x 600 x 380
12737839	-	TX150-ST5	5	Stainless steel	150 x 300 x 105	180 x 330 x 395
12747839	-	TX150-ST12	12	Stainless steel	300 x 325 x 150	330 x 360 x 395
12757839	-	TX150-ST18	18	Stainless steel	300 x 505 x 105	330 x 540 x 395
12767839	-	TX150-ST26	26	Stainless steel	300 x 505 x 200	330 x 540 x 395
12777839	-	TX150-ST38	38	Stainless steel	300 x 690 x 200	333 x 730 x 450

TXF200 thermostat, temperature range -50°C to 200°C

Catalogue No	Legacy No	Alt. No	Capacity, L	Bath material	Dimensions, internal [w x d x h], mm	Dimensions, external [w x d x h], mm
12747849	-	TXF200-P5	5	Plastic	160 x 240 x 150	200 x 390 x 380
12757849	-	TXF200-P12	12	Plastic	280 x 325 x 150	350 x 415 x 380
12767849	-	TXF200-P18	18	Plastic	290 x 510 x 150	350 x 600 x 380
12777849	-	TXF200-ST5	5	Stainless steel	150 x 300 x 105	180 x 330 x 395
12787849	-	TXF200-ST12	12	Stainless steel	300 x 325 x 150	330 x 360 x 395
12797849	-	TXF200-ST18	18	Stainless steel	300 x 505 x 105	330 x 540 x 395
12707859	-	TXF200-ST26	26	Stainless steel	300 x 505 x 200	330 x 540 x 395
12717859	-	TXF200-ST38	38	Stainless steel	300 x 690 x 200	333 x 730 x 450

Accessories

Catalogue No	Legacy No	Alt. No	Description
12787779	-	STL5	Flat lid, no hinges, stainless steel
12757779	-	STL12	Gabled lid, with removable hinges stainless steel
12767779	-	STL26	Gabled lid, with removable hinges stainless steel
12777779	-	STL38	Gabled lid, with removable hinges stainless steel
11436467	-	PL5	Curved lid, plastic
11446467	-	PL12	Curved lid, plastic
11456467	-	PL18	Curved lid, plastic

