

Baths, water, unstirred, accessories

Clifton Range

Test tube rack

Overall dimensions of rack 270mm x 70mm x 138mm

Maximum number of racks per bath

2L, 4L	1
8L	2
5L, 15L	4
9L, 22L, 28L	6
56L	12

Thermometer clip

Lies flush with top edge of bath to prevent accidental damage

Scale: 0 to 100°C in 2°C graduations

Lids

Catalogue No	Legacy No	Alt. No	Description
11419519	-	GL1-4	Gabled lid, hinged, for 4L baths
12758888	-	GL1-8	Gabled lid, hinged, clear, for 8L baths
12768888	-	GL1-14	Gabled lid, hinged, for 14L baths
6 11429519	-	SL1-4	Gabled lid, stainless steel, for 4L baths
11439519	-	SL1-8	Gabled lid, stainless steel, for 8L baths
11449519	-	SL1-14	Gabled lid, stainless steel, for 14L baths
11786577	-	GL1-22	Gabled lid, transparent
1 12778888	-	SL1-22	Gabled lid, stainless steel, for 9L, 22L or 28L baths
3 11894313	-	SL1-22H	Hinged gabled lid, stainless steel, for 22L or 28L baths
13121303	-	LD-4	Flat lid, for 4L baths
13131303	-	LD-8	Flat lid, for 8L baths
13141303	-	LD-14	Flat lid, for 14L baths
11766627	-	LD-22	Flat lid, for 9L, 22L, or 28L baths
11746597	-	SLR1-4	Stainless steel ring set 2 x 105mm holes for 4L baths
11429549	-	SLR1-8	Stainless steel ring set lid 4 x 83mm holes for 8L baths
2 11439549	-	SLR1-14	Stainless steel ring set lid 4 x 105mm holes for 14L baths
11449549	-	SLR2-14	Stainless steel ring set lid 6 x 83mm holes for 14L baths
12788888	-	SLR1-22	Stainless steel ring set lid 6 x 105mm holes for 9L, 22L or 28L baths

Polypropylene spheres

Catalogue No	Legacy No	Alt. No	Description	Pack qty
4 11469489	-	BP0368	Polypropylene spheres	200

Test tube racks

Catalogue No	Legacy No	Alt. No	Description
5 11419529	-	6870	Test tube rack, stainless steel, 26 x 17mm diameter
11439529	-	6871	Test tube rack, stainless steel, 16 x 26mm diameter
11409529	-	6872	Test tube rack, stainless steel, 36 x 13mm diameter
11429529	-	6873	Test tube rack, stainless steel, 18 x 19mm diameter
11449529	-	6900	Test tube rack, stainless steel, 12 x 32mm diameter

Catalogue No	Legacy No	Alt. No	Description
11409539	-	TC-1	Thermometer clip and thermometer
11419539	-	BX0616	Draining siphon
12778898	-	ST-1	Sterilising tray
12798898	-	SA01409	Pair of guide rail assemblies for sterilising tray



1



2



3



4



5



6