

Bottles

glass reagent



Bottles, premium, borosilicate glass, reagent and media, screw neck, ISO 4796 with retrace code



Ideal for demanding sterilisation processes such as hot air, dry sterilisation and depyrogenisation.

- Broader temperature resistant screw cap and matching pouring ring. The closure is fitted with a PTFE coated silicone seal and can withstand temperatures between -196°C to 260°C
- Leak free with regards to liquid and gaseous media
- No colourants in the closure or pouring ring, avoiding the risk of contamination to ovens, drying chambers and sterilisation cabinets, caused by the release of colourants
- Retrace code batch identification, which can be used to check details of manufacturing date, glass quality and conformity with standards. Batch certificate of traceability supplied with each pack
- Graduated scale accurate to $\pm 5\%$, plus extra subdivisions. Additional reverse scale permits the reading of current contents and quantities removed at the same time
- Guaranteed thermal shock resistance of 160K. The bottles have been tested in accordance with the strict guidelines of TUV Rheinland
- The whole system comprising bottle, closure and pouring ring conform to USP and FDA standards

Catalogue No	Legacy No	Alt. No	Capacity, mL	ISO thread GL	Height, mm	O.D., mm	Pack qty
13167644	-	1127075	100	45	105	56	10
13127444	-	1127076	250	45	143	70	10
13137444	-	1127077	500	45	181	86	10
13177644	-	1127078	1,000	45	230	101	10



Bottles, borosilicate glass, reagent and media screw neck, with retrace code



Range of clear and amber bottles with excellent chemical resistance, ISO thread, and graduations. Supplied without pouring ring and cap. Autoclavable.

For caps and pouring rings, see page 175.

Clear

Catalogue No	Legacy No	Alt. No	Capacity, mL	Height, mm	O.D., mm	ISO thread GL	Pack qty
13147344	-	2180114	25	70	36	25	10
13174784	-	2180117	50	87	46	32	10
13147154	-	2180124	100	100	56	45	10
13157246	-	218012906	150	110	62	45	10
13124324	-	2180136	250	138	70	45	10
13154574	-	2180144	500	176	86	45	10
13167246	-	218015106	750	203	95	45	10
13181744	-	2180154	1,000	225	105	45	10
12821114	-	2180163	2,000	260	136	45	1
13177246	-	218016908	3,500	295	160	45	1
13127264	-	2180173	5,000	330	186	45	1
13164294	-	2180186	10,000	410	234	45	1
13137584	-	2180188	15,000	445	268	45	1
13167913	-	2180191	20,000	505	300	45	1

Amber

Catalogue No	Legacy No	Alt. No	Capacity, mL	Height, mm	O.D., mm	ISO thread GL	Pack qty
13194594	-	2180614	25	70	36	25	10
13104344	-	2180617	50	87	46	32	10
13178023	-	2180624	100	100	56	45	10
13117246	-	218062902	150	110	62	45	10
13114344	-	2180636	250	138	70	45	10
13121764	-	2180644	500	176	86	45	10
13127246	-	218065102	750	203	95	45	10
12307198	-	2180654	1,000	225	105	45	1
13137246	-	218066304	2,000	260	136	45	1
13147246	-	218066904	3,500	295	160	45	1