

Bottles, polypropylene, aspirators, Thermo Scientific Nalgene®, Lowboy™



Excellent chemical resistance with leakproof closure and tap.

Autoclavable. Before autoclaving, closure should be disengaged from the threads.

Catalogue No	Legacy No	Alt. No	Capacity, L	Neck size, mm	Height, mm	O.D., mm	Closure size, mm
10799091	-	2324-0008	8	44	165	368	63
10662561	-	2324-0015	15	44	203	368	63

Spare part

Catalogue No	Legacy No	Alt. No	Description
10652561	-	6422-0010	Tap



Bottles, polypropylene, aspirators



Azlon

Squat form aspirators, for easy placement on shelves or bench top, with moulded in handle for safe and secure handling.

Moulded graduations in litres and moulded caps with drill guides for easy location of probes and feed lines. Stopcock opens and closes in a quarter turn. Autoclavable. Closures and stopcocks can be removed or loosened before autoclaving to avoid accidental implosion.

Catalogue No	Legacy No	Alt. No	Capacity, L	Cap size, mm	Height with cap, mm	Width, mm
11393454	-	BNP05A	5	45	270	220
11303464	-	BNP10A	10	90	360	280
11313464	-	BNP20A	20	90	420	350

Spare parts

Catalogue No	Legacy No	Alt. No	Description
11323464	-	CGZ85PP	Cap for 11303464, 11313464, 11383494, 11393494
11333464	-	SCBNP01	Spare stopcock



Bottles, polycarbonate, aspirators, Thermo Scientific Nalgene®, Clearboy™



Transparent and tough, with polypropylene tap and screw closure. Useful for many applications, including use in cold rooms of sub-zero temperatures.

Graduated in 5L increments. Supplied with tap closure. Leakproof. Autoclavable.

Catalogue No	Legacy No	Alt. No	Capacity, L	Neck size, mm	Dimensions [O.D. x h], mm	Closure size, mm
10722332	-	2317-0020	10	64	249 x 389	83B
10684071	-	2317-0050	20	64	282 x 513	83B

Spare part

Catalogue No	Legacy No	Alt. No	Description
10652561	-	6422-0010	Tap

