



Bottles, HDPE, translucent, Thermo Scientific Nalgene®



- Bottles and closures are supplied assembled to minimise contamination during transit and storage
- Leakproof
- Excellent resistance to most acids, bases and alcohols
- Narrow neck is ideal for pouring liquids
- Wide neck is easy to fill with dry materials or liquids
- HDPE with polypropylene screw closure

Break resistant bottles in narrow and wide neck configurations for a variety of field sampling needs.

Narrow neck

Catalogue No	Legacy No	Alt. No	Capacity, mL	Neck size, mm	Height, mm	O.D., mm	Closure size, mm	Pack qty
10441781	-	2089-0001	30	14	58	36	20	72
10490333	-	2089-0002	60	14	81	38	20	72
10253971	-	2089-0004	125	17	97	48	24	72
10400343	-	2089-0008	250	17	130	61	24	72
10601722	-	2089-0016	500	21	171	71	28	48
10653252	-	2089-0032	1,000	27	213	91	38	24

Narrow neck, bulk pack

Catalogue No	Legacy No	Alt. No	Capacity, mL	Neck size, mm	Height, mm	O.D., mm	Closure size, mm	Pack qty
10334912	-	33-2089-0004	125	17	97	48	24	500
10035260	-	33-2089-0008	250	17	130	61	24	250
10209090	-	33-2089-0012	500	21	171	71	28	125
10294351	-	33-2089-0032	1,000	27	213	91	38	50

Wide neck

Catalogue No	Legacy No	Alt. No	Capacity, mL	Neck size, mm	Height, mm	O.D., mm	Closure size, mm	Pack qty
10214071	-	2189-0001	30	20	64	36	28	72
10654212	-	2189-0002	60	20	86	38	28	72
10102881	-	2189-0004	125	28	99	51	38	72
10664212	-	2189-0008	250	33	132	62	43	72
10244351	-	2189-0016	500	43	168	74	53	48
10505021	-	2189-0032	1,000	51	201	91	63	24

Wide neck, bulk pack

Catalogue No	Legacy No	Alt. No	Capacity, mL	Neck size, mm	Height, mm	O.D., mm	Closure size, mm	Pack qty
10219120	-	33-2189-0002	60	20	86	38	28	1,000
10401571	-	33-2189-0004	125	28	99	51	38	500
10300441	-	33-2189-0008	250	33	132	62	43	250
10267780	-	33-2189-0016	500	43	168	74	53	125