

# Bottles HDPE



## Bottles, HDPE, translucent, narrow neck, Thermo Scientific Nalgene®



- Good freezer storage to -100°C
- All purpose, highly reliable and durable
- For long term use
- Narrow neck, ideal for storing samples
- Can be used with most corrosives
- Amber bottles help protect and prolong life of light-sensitive materials

Bottles with excellent chemical resistance. Ideal for storing, shipping and packaging liquids.

### Clear

| Catalogue No | Legacy No | Alt. No   | Capacity, mL | Neck size, mm | Height, mm | O.D., mm | Closure size, mm | Pack qty |
|--------------|-----------|-----------|--------------|---------------|------------|----------|------------------|----------|
| 10718331     | -         | 2002-9125 | 4            | 9             | 41         | 16       | 13               | 12       |
| 10550022     | -         | 2002-9025 | 8            | 13            | 45         | 23       | 20               | 12       |
| 10758711     | -         | 2002-9050 | 15           | 13            | 56         | 23       | 20               | 12       |
| 10768711     | -         | 2002-0001 | 30           | 13            | 61         | 36       | 20               | 12       |
| 10313501     | -         | 2002-0002 | 60           | 13            | 84         | 38       | 20               | 12       |
| 10159140     | -         | 2002-0004 | 125          | 18            | 102        | 51       | 24               | 12       |
| 10782711     | -         | 2002-0006 | 175          | 18            | 125        | 54       | 24               | 12       |
| 10507341     | -         | 2002-0008 | 250          | 18            | 135        | 61       | 24               | 12       |
| 10517341     | -         | 2002-0016 | 500          | 23            | 170        | 74       | 28               | 12       |
| 10027950     | -         | 2002-0032 | 1,000        | 28            | 216        | 91       | 38               | 6        |

### Amber

| Catalogue No | Legacy No | Alt. No   | Capacity, mL | Neck size, mm | Height, mm | O.D., mm | Closure size, mm | Pack qty |
|--------------|-----------|-----------|--------------|---------------|------------|----------|------------------|----------|
| 10119960     | -         | 2004-9125 | 4            | 9             | 41         | 16       | 13               | 12       |
| 10142641     | -         | 2004-9025 | 8            | 13            | 45         | 23       | 20               | 12       |
| 10749481     | -         | 2004-9050 | 15           | 13            | 56         | 23       | 20               | 12       |
| 10474481     | -         | 2004-0001 | 30           | 13            | 61         | 36       | 20               | 12       |
| 10343691     | -         | 2004-0002 | 60           | 13            | 84         | 38       | 20               | 12       |
| 10588101     | -         | 2004-0004 | 125          | 18            | 102        | 51       | 24               | 12       |
| 10140101     | -         | 2004-0008 | 250          | 18            | 135        | 61       | 24               | 12       |
| 10281941     | -         | 2004-0016 | 500          | 23            | 170        | 74       | 28               | 12       |
| 10057960     | -         | 2004-0032 | 1,000        | 28            | 216        | 91       | 38               | 6        |



## Bottles, HDPE, translucent, medium weight, narrow neck



Ideal for general liquid storage and sampling.

Supplied capped and in shrink-wrapped packs. All sizes feature S.P.I. recycling codes on their base. Polypropylene screw closure.

| Catalogue No | Legacy No | Alt. No   | Capacity, mL | Neck size, mm | Height, mm | O.D., mm | Pack qty |
|--------------|-----------|-----------|--------------|---------------|------------|----------|----------|
| 10608015     | -         | BNH0030PN | 30           | 17            | 69         | 32       | 10       |
| 10195714     | -         | BNH0060PN | 60           | 17            | 81         | 40       | 10       |
| 10616865     | -         | BNH0125PN | 125          | 22            | 104        | 50       | 10       |
| 10626865     | -         | BNH0150N  | 150          | 22            | 116        | 50       | 10       |
| 12074707     | -         | BNH0250P  | 250          | 22            | 140        | 60       | 10       |
| 10708805     | -         | BNH0500P  | 500          | 22            | 161        | 75       | 10       |
| 10125374     | -         | BNH1000P  | 1000         | 22            | 201        | 90       | 5        |
| 10298144     | -         | BNH2000P  | 2,000        | 37            | 242        | 120      | 5        |