

Bottles HDPE



Bottles, HDPE, amber, wide neck, Fisherbrand®



Fisherbrand®

- Opaque amber reduces UV light transmissions to protect light sensitive contents
- Wide neck allows easy retrieval of grab samples and reduces turbulence when filling
- Excellent chemical resistance to most acids, bases and alcohol
- Leakproof

Catalogue No	Legacy No	Capacity, mL	Pack qty
12662415	-	60	72
12672415	-	125	72
12682415	-	250	72
12602425	-	500	48
12692415	-	1,000	24



Bottles, HDPE, translucent, wide neck, Thermo Scientific Nalgene®

Thermo
SCIENTIFIC | **NALGENE**

- Suitable for shipping liquids
- Polypropylene screw closure
- Leakproof
- Excellent durability and sample protection
- More rigid than LDPE
- Wide mouth, easy to fill with dry materials or liquids
- Excellent chemical resistance to most acids, bases and alcohols
- Good for freezer use up to -100°C
- Amber bottles help protect and prolong life of light-sensitive materials

Clear and amber durable, general purpose bottles.



Clear

Catalogue No	Legacy No	Alt. No	Capacity, mL	Neck I.D., mm	Height, mm	O.D., mm	Closure size, mm	Pack qty
10353501	-	2104-0001	30	23	63	36	28	12
10559831	-	2104-0002	60	20	86	38	28	12
10569831	-	2104-0004	125	31	99	51	38	12
10623141	-	2104-0008	250	33	131	61	43	12
10633141	-	2104-0016	500	43	168	74	53	12
10547341	-	2104-0032	1,000	53	199	91	63	6

Amber

Catalogue No	Legacy No	Alt. No	Capacity, mL	Neck I.D., mm	Height, mm	O.D., mm	Closure size, mm	Pack qty
10547692	-	2106-0002	60	20	86	39	28	12
10024790	-	2106-0004	125	31	99	50	38	12
10034790	-	2106-0008	250	33	131	61	43	12
10665601	-	2106-0016	500	43	168	73	53	12
10064810	-	2106-0032	1,000	53	199	91	63	6