

Bottles, polypropylene, translucent, Fisherbrand®



Fisherbrand

- Translucent, leakproof
- Excellent chemical resistance
- HDPE closure
- Autoclavable

Catalogue No	Legacy No	Capacity, mL	Neck size, mm	Height, mm	O.D., mm	Pack qty
11561933	-	250	45	120	58	10
11571933	-	500	45	150	75	10



Bottles, polypropylene, translucent, wide neck, Fisherbrand®



Fisherbrand

Fisherbrand® translucent bottles offer flexibility for countless applications.

- General purpose - fill with dry materials or liquids
- Chemically resistant to a wide range of acids and organic alcohols
- Wide neck allows easy retrieval of grab samples and reduces turbulence when filling

Catalogue No	Legacy No	Capacity, mL	Pack qty
12632405	-	30	72
12642405	-	60	72
12662405	-	125	72
12672405	-	250	72
12682405	-	500	48
12612415	-	1,000	24



Bottles, polypropylene (PPCO), translucent, narrow neck, Thermo Scientific Nalgene®



Thermo Scientific NALGENE

- Polypropylene screw closure
- Autoclavable with cap loosened off threads
- Excellent for chemicals, specimens or materials which require autoclaving
- Withstands temperatures up to +121°C
- PPCO (polypropylene co-polymer) offers better clarity than HDPE or LDPE
- Narrow neck ideal for storing, shipping and packing liquids
- Resistant to most acids, bases and alcohols

Ideal for general laboratory applications.

Catalogue No	Legacy No	Capacity, mL	Neck size, mm	Height, mm	O.D., mm	Closure size, mm	Pack qty
10150391	-	4	8	40	16	13	12
10444681	-	8	14	45	25	20	12
10241901	-	15	13	56	25	20	12
10779481	-	30	13	61	36	20	12
10363691	-	60	13	84	38	20	12
10190291	-	125	18	102	51	24	12
10159760	-	250	18	132	61	24	12
10443151	-	500	23	170	74	28	12
10363401	-	1,000	28	216	91	38	6

