

Bottles polypropylene

Bottles, polypropylene (PPCO), economy, Thermo Scientific Nalgene®



Suitable for short term general purpose laboratory applications.

- Made from PPCO (polypropylene co-polymer)
- Lightweight
- Narrow neck for easy pouring
- Wide neck is easy to fill with liquids or dry materials
- Excellent chemical resistance to most acids, bases and alcohols
- Can be autoclaved

Narrow neck

Catalogue No	Legacy No	Alt. No	Capacity	Neck size, mm	Height, mm	O.D., mm	Closure size, mm	Pack qty
11399604	-	2087-0001	30	14	58	36	20	72
11309614	-	2087-0002	60	14	81	38	20	72
11319614	-	2087-0004	125	17	97	48	24	72
11329614	-	2087-0008	250	17	130	61	24	72
11339614	-	2087-0016	500	21	171	71	28	48
11349614	-	2087-0032	1,000	21	213	91	38	24

Wide neck

Catalogue No	Legacy No	Alt. No	Capacity, mL	Neck size, mm	Height, mm	O.D., mm	Closure size, mm	Pack qty
11399614	-	2187-0001	30	20	64	36	28	72
11309624	-	2187-0002	60	20	86	38	28	72
11319624	-	2187-0004	125	28	99	51	38	72
11329624	-	2187-0008	250	33	132	62	43	72
11339624	-	2187-0016	500	43	168	74	53	48
11349624	-	2187-0032	1,000	51	201	91	63	24



Bottles, polypropylene, translucent, media, narrow neck



Azlon

Ideal for liquid storage, media preparation, etc.

These high quality bottles are supplied capped, in shrinkwrapped packs. With S.P.I. recycling codes. Autoclavable. Excellent general chemical resistance.

Catalogue No	Legacy No	Alt. No	Capacity, mL	Height, mm	O.D., mm	Neck size, mm	Pack qty
10044514	-	BNP0030PN	30	69	32	17	10
10004564	-	BNP0060PN	60	81	40	17	10
10278304	-	BNP0125PN	125	104	50	22	10
10472395	-	BNP0150PN	150	116	50	22	10
10175564	-	BNP0250P	250	140	60	22	10
10490655	-	BNP0500P	500	161	75	22	10
10207774	-	BNP1000P	1,000	201	90	22	5
10419704	-	BNP2000P	2,000	242	120	37	5