

Centrifuge tubes



Tubes, centrifuge, conical, 15mL/50mL, Fisherbrand®



- A large white marking area makes it easy to label samples
- All sizes and types are available with a blue plug seal closure with a double start thread design that provides a tight secure seal
- Most types are also available with a standard flat top cap
- Tubes are sterilised by gamma radiation and are non-cytotoxic
- Available packed in racks or sleeves
- Racks are not available separately

Disposable conical bottom tubes with evenly spaced black printed graduations to enhance readability without losing precision.

15mL tubes, graduated in 0.5mL subdivisions from 2mL to 14.5mL, 17mm x 119mm

Polypropylene, packed in racks, sterile

Catalogue No	Legacy No	Closure style	RCF, x g [max.]	Inner pack qty	Pack qty
11542293	-	Plug seal	6,000	50	500
11849650	-	Flat top	6,000	50	500

Polyethylene terephthalate (PET), packed in racks

Catalogue No	Legacy No	Closure style	RCF, x g [max.]	Inner pack qty	Pack qty
11879640	-	Plug seal	1,800	50	500

50mL tubes, graduated in 5mL subdivisions from 5mL to 50mL, 28mm x 115mm

Polypropylene, all are sterile except 11829650

Catalogue No	Legacy No	Closure style	RCF, x g [max.]	Pack type	Inner pack qty	Pack qty
11899640	-	Plug seal	6,000	Rack	25	500
11809650	-	Plug seal	6,000	Sleeve	25	500
11512303	-	Flat top	6,000	Sleeve	25	500
11819650	-	Flat top	6,000	Rack	25	500
11829650	-	Plug seal	6,000	Sleeve	25	500

Polyethylene terephthalate (PET), packed in racks

Catalogue No	Legacy No	Closure style	RCF, x g [max.]	Inner pack qty	Pack qty
11839650	-	Plug seal	1,800	25	500



Tubes, centrifuge, Thermo Scientific Nunc® EZ Flip™ cap



- Integrated leakproof cap prevents cross contamination and lost closures
- Patented hinged cap system
- Ergonomic design for one handed opening and closing
- Clarified polypropylene construction provides chemical resistance and high RCF tolerance
- White writing area and flat caps for sample identification
- Printed graduations
- Standard dimensions compatible with centrifuge rotors
- 15mL and 50mL versions are available

Thermo Scientific Nunc EZ Flip™ conical centrifuge tubes feature a patented, hinged flip cap design in standard 15mL and 50mL sizes.

These innovative tubes can be opened and closed with one hand, freeing the user's other hand for manipulations, such as pipetting.

Catalogue No	Legacy No	Alt. No	Description	RCF, x g [max.]	Inner pack qty	Pack qty
10512802	-	362694	15mL bulk	8,500	50	500
10777414	-	362695	15mL racked	8,500	50	500
10048631	-	362696	50mL bulk	9,500	25	500
10571663	-	362697	50mL racked	9,500	25	500