

Vials, crimp, 20mm opening, headspace, Fisherbrand®


Fisherbrand®



Catalogue No	Legacy No	Description	For use with	Capacity, mL	Dimensions, mm	Pack qty
11510545	-	Headspace vial, crimp, clear glass, flat bottom	Varian	5	20 x 38	100
12931231	-	Headspace vial, crimp, clear glass, round bottom	Perkin Elmer	5	22 x 38.2	100
11520545	-	Headspace vial, crimp, clear glass, flat bottom	Varian	10	20 x 54.5	100
12941231	-	Headspace vial, DIN-crimp, clear glass, round bottom	Carlo Erba, CTC, Fisons, Varian (CP)	10	22.5 x 46	100
12951231	-	Headspace vial, DIN-crimp, clear glass, long neck, flat bottom	Carlo Erba, Dani, Fisons, Agilent	10	22.5 x 46	100
11540535	-	Headspace vial, crimp, clear glass, round bottom	Perkin Elmer, Tekmar	20	23 x 75.5	100
12961231	-	Headspace vial, crimp, clear glass, flat bottom, graduated, with marking spot	Perkin Elmer, Tekmar	20	23 x 75.5	100
12971231	-	Headspace vial, DIN-crimp, clear glass, long neck, flat bottom	Carlo Erba, Dani, Fisons, Agilent	20	22.5 x 75.5	100
12981231	-	Headspace vial, DIN-crimp, clear glass, long neck, round bottom	CTC PAL (Varian, Gerstel, Atas, Shimadzu) + TriPlus HS)	20	22.5 x 75.5	100
11530535	-	Headspace vial, crimp, clear glass, flat bottom	-	50	31 x 101	100
12921231	-	Headspace vial, DIN-crimp, amber glass, round bottom	Carlo Erba, CTC, Fisons, Varian (CP)	10	22.5 x 46	100
12981241	-	Headspace vial, crimp, amber glass, round bottom	Perkin Elmer	20	23 x 75.5	100
12910991	-	Headspace vial, DIN-crimp, amber glass, long neck, round bottom	CTC PAL (Varian, Gerstel, Atas, Shimadzu) + TriPlus HS)	20	22.5 x 75.5	100

Caps, crimp, aluminium, 20mm, headspace, Fisherbrand®

Fisherbrand®



Catalogue No	Legacy No	Description	Septa material	Thickness, mm	Pack qty
12991231	-	Aluminium cap, plain, 10mm centre hole	Butyl septa, non-PTFE laminated	3.0	100
12901241	-	Headspace cap, clear, lacquered	Butyl septa, non-PTFE laminated	3.0	100
11510515	-	Magnetic cap, gold, 8mm centre hole	Butyl septa, non-PTFE laminated	3.0	100
10132712	-	Aluminium cap, plain, 10mm centre hole	Butyl/grey PTFE septa, all PTFE laminated	3.0	100
12921241	-	Headspace cap, clear lacquered	Butyl/grey PTFE septa, all PTFE laminated	3.0	100
11520515	-	Magnetic cap, gold, 8mm centre hole	Butyl/grey PTFE septa, all PTFE laminated	3.0	100
12931241	-	Aluminium cap, plain, 10mm centre hole	With Pharma-Fix butyl/PTFE septa	3.0	100
12941241	-	Headspace cap, clear, lacquered	With Pharma-Fix butyl/PTFE septa	3.0	100
11530515	-	Magnetic cap, gold, 8mm centre hole	With Pharma-Fix butyl/PTFE septa	3.0	100
11506144	-	Aluminium cap, plain, 10mm centre hole	Transparent blue silicone/white PTFE, all PTFE laminated	3.0	100
12951241	-	Headspace cap, clear, lacquered	Transparent blue silicone/white PTFE, all PTFE laminated	3.0	100
12961241	-	Magnetic cap, gold, 8mm centre hole	Transparent blue silicone/transparent PTFE, all PTFE laminated	3.0	100
10617445	-	Magnetic bimetal cap, red, 8mm centre hole	Transparent blue silicone/transparent PTFE, all PTFE laminated	3.0	100
11590505	-	Aluminium cap, plain, 10mm centre hole	White silicone/beige PTFE, all PTFE laminated	3.25	100
11580505	-	Headspace cap, clear, lacquered	White silicone/beige PTFE (HT Agilent quality), all PTFE laminated	3.25	100
12990971	-	Aluminium cap, plain, centre hole	Silicone white/aluminium foil silver	3.0	100
12980971	-	Headspace cap, clear lacquered, centre hole	Silicone white/aluminium foil silver	3.0	100
12930981	-	Magnetic cap, gold, 5mm centre hole	Silicone white/aluminium foil silver	3.0	100
12970981	-	Magnetic Bi-Matall cap, red, 8mm centre hole	Butyl/PTFE grey	3.0	100
11368641	-	Magnetic Bi-Matall cap, red, 8mm centre hole	Silicone white/PTFE beige	3.2	100
11540525	-	Magnetic cap, gold, 8mm centre hole	Silicone white/PTFE beige	3.2	100