





Storage boxes, cryogenic, polypropylene, with grid inserts, Fisherbrand®



Fisherbrand®

For sample storage and transportation.


- Autoclavable, temperature resistant from -90°C to +121°C
- High quality polypropylene
- Available in natural and various colours
- 81 places, 9 x 9 grid, numerically coded
- Robust hinge with safe snap-on lid
- Stackable for safe transportation
- Drain bores at the bottom for dew liquids
- Removable grid inserts
- Can also be used for larger vessels by removing the grid inserts

Dimensions (l x w x h), 133mm x 133mm x 52mm, 81 place, 9 x 9 grid

Catalogue No	Legacy No	Accommodates	Colour	Pack qty
11587713	-	81 x 1.2mL to 2.0mL vials	Natural	5
11597713	-	81 x 1.2mL to 2.0mL vials	Yellow	5
11507723	-	81 x 1.2mL to 2.0mL vials	Red	5
11517723	-	81 x 1.2mL to 2.0mL vials	Green	5
11527723	-	81 x 1.2mL to 2.0mL vials	Blue	5

Dimensions (l x w x h), 133mm x 133mm x 75mm, 81 place, 9 x 9 grid

Catalogue No	Legacy No	Accommodates	Colour	Pack qty
11866903	-	81 x 3.0mL to 5.0mL vials	Natural	5
11876903	-	81 x 3.0mL to 4.0mL vials	Yellow	5
11886903	-	81 x 3.0mL to 5.0mL vials	Red	5
11896903	-	81 x 3.0mL to 5.0mL vials	Green	5
11806913	-	81 x 3.0mL to 5.0mL vials	Blue	5



Storage boxes, cryogenic, polycarbonate, Fisherbrand®

Fisherbrand®

- Polycarbonate
- Separate lid
- Grid reference printed on lid and base
- Can be used in mechanical freezers and liquid nitrogen dewars

Storage boxes, printed lid, 81 place, 9 x 9 grid

Catalogue No	Legacy No	Description	Dimensions [l x w x h], mm	Pack qty
12902431	-	Storage box for 1.2mL to 2.0mL vials	133 x 133 x 52	24
12912431	-	Storage box for 5.0mL vials	133 x 133 x 95	24
12932431	-	Storage box for 1.2mL to 2.0mL vials	133 x 133 x 52	4
12922431	-	Storage box for 5.0mL vials	133 x 133 x 95	4