



Storage boxes, cryogenic, Thermo Scientific Nalgene® CryoBox™



For ultra-low temperature storage of cryogenic vials and other vials; ideal for use in mechanical and liquid nitrogen freezers.

- Made of polycarbonate, a durable and economic alternative to cardboard or stainless steel
- Usable temperature range of -196°C to 121°C
- Feature an interior grid with moulded numbers
- Lid is printed in blue (except for 10038700, in black) with a grid of corresponding numbers
- Lids keyed to one position so grids always match
- Vented to allow rapid thawing
- Lid can easily be removed straight from the freezer

Catalogue No	Legacy No	Alt. No	Array	Dimensions [w x d x h], mm	Capacity, mL	Colour
10739331	-	5025-0505	5 x 5	76 x 76 x 51	2 (vial)	White
10324741	-	5026-0909	9 x 9	133 x 133 x 51	2 (vial)	White
10038700	-	5027-0909	9 x 9	133 x 133 x 95	2 (vial)	White
11754703	-	867013-0240	9 x 9	133 x 133 x 51	2 (vial)	Red grid
11764703	-	867013-0241	9 x 9	133 x 133 x 51	2 (vial)	Yellow grid
11391514	-	867013-0242	9 x 9	133 x 133 x 51	2 (vial)	Green grid
11301524	-	867013-0243	9 x 9	133 x 133 x 51	2 (vial)	Blue grid
11311524	-	867013-0244	9 x 9	133 x 133 x 51	2 (vial)	Grey grid



Storage racks, freezer, Fisherbrand®



- Corrosion resistant stainless steel
- Accommodate a variety of racks and storage boxes
- Side opening with handles on the front for easy removal
- Designed to fit a wide range of freezer brands and makes

	Catalogue No	Legacy No	Description	Dimensions [w x d x h], mm	For use with
1	12653476	-	Holds 4 x 11700634 15mL and 50mL storage boxes	154 x 157 x 539 (external) 152 x 156 x 133 (each slot)	Chest freezers
2	12643476	-	Holds 6 x 11700634 15mL and 50mL storage boxes	459 x 159 x 271 (external) 152 x 157 x 134 (each slot)	Upright freezers
3	12902690	-	Holds 16 x PCR tube racks (11710344), 48 x 96 well assay plates, 32 x 24 well assay plates, 16 deep well plates	543 x 112 x 330 (external) 135 x 111 x 56 (each slot)	Upright freezers
4	12912690	-	Holds 18 x microscope slide boxes (208mm x 175mm x 34mm), 9 x reversible PCR racks (11787265), 9 x 50 well EPS tube racks (11948084)	541 x 217 x 221 (external) 180 x 212 x 73 (each slot)	Upright freezers

