


Cylinders

borosilicate glass Class A



Cylinders, borosilicate glass, spouted, Class A, blue ceramic graduations, Fisherbrand®

Fisherbrand®

DIN 12680 BS 604 ISO 4788

Graduated in blue ceramic. Hexagonal base. Individual reorder code on each item.

Catalogue No	Legacy No	Capacity, mL	Graduations, mL	Pack qty
12952310	-	5	0.1	2
11517832	-	10	0.2	2
11527832	-	25	0.5	2
11537832	-	50	1	2
11547832	-	100	1	2
11557832	-	250	2	2
11577832	-	1,000	10	2
12962320	-	2,000	20	2



Cylinders, borosilicate glass, spouted, Class A, white enamel graduations

PYREX®

DIN 12680, BS 604 & ISO 4788

Manufactured from Pyrex® borosilicate glass with graduations in white enamel for improved visibility with dark colours. Hexagonal base. Available with or without a certificate of calibration. Certified cylinders have an individual serial number for traceability to national standards, and a certificate showing the exact volume contained at 5 points on the scale. Batch certificate is downloadable.

Without certificate

Catalogue No	Legacy No	Alt. No	Capacity, mL	Graduations, mL	Tolerance, mL	Pack qty
11903789	-	3201/02	5	0.1	0.05	2
12594277	-	3201/04	10	0.2	0.10	2
11913789	-	3201/06	25	0.5	0.25	2
11923789	-	3201/08	50	1.0	0.50	2
11933789	-	3201/10	100	1.0	0.50	2
12093333	-	3201/12	250	2.0	1.0	2
11943789	-	3201/14	500	5.0	2.5	2
10793252	-	3201/16	1,000	10.0	5.0	1
12049867	-	3201/18	2,000	20.0	10.0	1

With certificate of calibration

Catalogue No	Legacy No	Alt. No	Capacity, mL	Graduations, mL	Tolerance, mL	Pack qty
11383544	-	3201/02WC	5	0.1	0.05	2
11393544	-	3201/04WC	10	0.2	0.10	2
11303554	-	3201/06WC	25	0.5	0.25	2
11313554	-	3201/08WC	50	1.0	0.50	2
11323554	-	3201/10WC	100	1.0	0.50	2
11333554	-	3201/12WC	250	2.0	1.0	2
11343554	-	3201/14WC	500	5.0	2.5	2
11353554	-	3201/16WC	1,000	10.0	5.0	1
11363554	-	3201/18WC	2,000	20.0	10.0	1

