

Cylinders

polypropylene



Cylinders, polypropylene, translucent, spouted, moulded graduations, Thermo Scientific Nalgene®



Moulded graduations, non-wetting and large base resists tipping. Hexagonal base. Can be chemically sterilised, which will not affect accuracy. Autoclavable, will have some effect on accuracy.

Catalogue No	Legacy No	Alt. No	Capacity, mL	Graduations, mL
10190181	-	3662-0010	10	0.2
10403921	-	3662-0025	25	0.5
10559101	-	3662-0050	50	1.0
10683931	-	3662-0100	100	1.0
10078510	-	3662-0250	250	2.0
10624121	-	3662-0500	500	5.0
10453921	-	3662-1000	1,000	10.0
10404091	-	3662-2000	2,000	20.0
10599481	-	3662-4000	4,000	50.0



Cylinders, polypropylene, translucent, spouted, blue graduations



ISO 6706.

Ultra clear, blue, printed graduations. Hexagonal base. With S.P.I. recycling code. Food safe.

Catalogue No	Legacy No	Alt. No	Capacity, mL	Graduations, mL	Pack qty
11373584	-	CPB0010P	10	0.2	10
11383584	-	CPB0025P	25	0.5	10
11393584	-	CPB0050P	50	1.0	10
11303594	-	CPB0100P	100	1.0	5
11940985	-	CPB0250P	250	2.0	1
11950985	-	CPB0500P	500	5.0	1
11333594	-	CPB1000P	1,000	10.0	1
11343594	-	CPB2000P	2,000	20.0	1



Cylinders, polypropylene, translucent, spouted, moulded graduations



ISO 6706.

Moulded graduations. Hexagonal base. With S.P.I. recycling code.

Catalogue No	Legacy No	Alt. No	Capacity, mL	Graduations, mL	Pack qty
11373574	-	CP0010P	10	0.2	10
11383574	-	CP0025P	25	0.5	10
10667045	-	CP0050P	50	1.0	10
10155234	-	CP0100P	100	1.0	5
11920985	-	CP0250P	250	2.0	1
11930985	-	CP0500P	500	5.0	1
10195184	-	CP1000P	1,000	10.0	1
11313584	-	CP2000P	2,000	20.0	1