

Extraction thimbles, glass microfibre, quartz microfibre

Whatman
Part of GE Healthcare



- Smoke stack gas monitoring
- Soxhlet extraction
- Analysing pesticide residues
- Determining oil/fat content of foods (e.g. French fries)
- Analysis of oil and grease in solid wastes
- Available in a range of sizes and wall thicknesses to suit your application
- Designed to fit most commercially available Soxhlet extractors

High-purity glass microfibre thimbles

The thimbles are completely free of binders or additives and can be used at temperatures up to 500°C or when using solvents incompatible with cellulose thimbles.

These thimbles are also used in pollution monitoring techniques (0.8µm nominal particle retention). Typical thickness 1.7 mm.

Quartz microfibre thimbles

Made from high purity quartz microfibre, this thimble is able to withstand high temperatures (up to 1,000°C). Suitable for both solvent extraction and air sampling applications.

Standard glass microfibre extraction thimbles

Grade 603 g (Borosilicate glass fibre with inorganic binder)

Catalogue No	Legacy No	Alt. No	Size, mm	Thickness, mm	Pack qty
10151141	-	10371005	16 x 50	1.0	25
10053431	-	10371007	19 x 90	1.0	25
10627491	-	10371011	22 x 80	1.5	25
10355131	-	10371019	25 x 100	1.5	25
10482794	-	10371023	26 x 100	1.5	25
10315181	-	10371025	28 x 60	1.5	25
11875335	-	10371029	25 x 98	1.5	25
10504602	-	10371036	30 x 100	1.5	25
10226981	-	10371042	33 x 94	1.5	25
10111241	-	10371043	33 x 100	1.5	25
10121241	-	10371045	33 x 118	1.5	25
10671852	-	10371055	35 x 150	1.5	25
11890775	-	10371103	10 x 38	1.5	25
10171331	-	10371114	26.8 x 68	1.5	25

High-purity glass microfibre extraction thimbles

Grade HP-GF (100% pure borosilicate glass)

Catalogue No	Legacy No	Alt. No	Size, mm	Pack qty
13178543	-	2814-199	19 x 90	25
11890755	-	2814-259	Tapered, 25 x 90	25
11894265	-	2814-300	30 x 100	25
13111914	-	2814-432	43 x 123	25
11874764	-	2814-533	53 x 135	25

Quartz microfibre (silica) extraction thimbles

Catalogue No	Legacy No	Alt. No	Size, mm	Pack qty
11840705	-	2812-259	Tapered, 25 x 90	10