

Incubator, shaking, bench top, ES-20

Grant bio

- Compact, bench top unit
- Simple digital programming of all parameters
- Clear indication of set and actual readings
- Full range of accessory platforms

A versatile and programmable bench top shaker-incubator for mixing and incubating biological cultures and samples in a variety of flasks and vessels. Equipped with direct drive shaking system for reliable, long term operation and forced air circulation ensuring constant and even temperature within the chamber.



Technical Specification

Load capacity, kg 2.5

Speed, rpm 50 to 250 (variable)

Temperature range, °C 25 to 42

Catalogue No	Legacy No	Alt. No	Description
11411058	-	ES-20	Shaker-incubator, bench top

Accessories

Catalogue No	Legacy No	Alt. No	Description
11461058	-	P6-250	Platform with clamps for 6 x 250mL flasks/400mL beakers
11451058	-	P12-100	Platform with clamps for 12 x 100mL flasks/150mL beakers
11441058	-	PP-4	Flat platform, with non-slip mat (219mm x 219 mm)

Choose online ordering

- Browse all our products and services with only a few clicks
- Take advantage of promotions throughout the year
- Orders are instantly and automatically processed
- Keep track of your budget using the order history functionality
- Access your contract prices
- Track your order online

UK and International: www.fisher.co.uk | **Ireland:** www.ie.fishersci.com

Sweden: www.fishersci.se | **Norway:** www.fishersci.no

Finland: www.fishersci.fi | **Denmark:** www.fishersci.dk

