

Inoculating loops, with needle, polystyrene, sterile, Fisherbrand®



Fisherbrand®

Sterilised at 10KGy, these inoculating needles/loops facilitate uniform and smooth streaking without damaging the gel surface.

With a loop at one end in either 1µL (crystal clear) or 10µL (blue) and a needle point at the other end, they are both available in three pack quantities depending upon the quantity in the inner pack. Shelf life printed on each sachet.

Catalogue No	Legacy No	Volume, µL	Colour	Inner pack qty	Pack qty
12662275	-	1	Clear	1	900
12682275	-	1	Clear	20	9,000
11753124	-	10	Blue	1	900
12612285	-	10	Blue	20	9,000
11773124	-	10	Blue	40	18,000



For pH indicator papers, see pages 605 to 607



Inoculating loops, polystyrene, sterile



High impact, polystyrene, disposable inoculating loops. 200mm long to reduce the risk of contamination.

Catalogue No	Legacy No	Alt. No	Capacity, µL	Inner pack qty	Pack qty
12870155	-	PL2010/1000	10	20	1,000
11576863	-	PL205/1000	5	20	1,000
12890155	-	PL201/1000	1	20	1,000



Inoculating loops, with sphere, sterile, Thermo Scientific Sterilin® Quadloops®




Thermo SCIENTIFIC

Sterile disposable inoculating device for either sampling fixed volumes or for dilution streaking to obtain isolated colonies.

Unique sphere/handle form provides up to four sterile streaking surfaces simply by turning the handle 90° after each streaking, thereby exposing a new sterile streaking surface.

Catalogue No	Legacy No	Alt. No	Capacity, µL	Colour	Pack qty
11778858	-	QL1	1	Green	1,000
11788858	-	QL10	10	Blue	1,000



Inoculating needles, polystyrene, sterile



Technical Service Consultants Ltd

Clear polystyrene with chisel tip to facilitate stab cultures and recovery of adherent colonies.

Catalogue No	Legacy No	Description	Inner pack qty	Pack qty
12391309	-	Clear polystyrene with chisel tip	20	1,000
12719885	-	Clear polystyrene with chisel tip individually wrapped	1	1,000

