


Microscope slide boxes, polystyrene

Kartell

Slots are well spaced for easy handling of 76mm x 26mm slides. Each slot is clearly numbered and can be cross referenced with corresponding index card. Lids are totally removable. Top and back clearly indicated. The 50 and 100 place boxes have locating lugs so boxes can be stacked. Made of virtually unbreakable high impact polystyrene.

Catalogue No	Legacy No	Capacity, slides	Width, mm	Depth, mm	Height, mm
11372874	-	25	98	83	38
10576191	-	50	230	97	35
10586191	-	100	230	180	35



Microscope slide staining tray, Fisherbrand®



Fisherbrand®

Stain, rinse and dry slides on a single working tray.

- Disposable microscope slide staining trays and lid (pack includes 4 base trays and 1 lid)
- Each tray holds up to eight slides and 38mL of liquid
- Compact size and recessed handles for easy transport
- Dark lid protects slides for light sensitive applications

Catalogue No	Legacy No	Material	Colour	Dimensions [w x d x h], mm
11968004	-	Polypropylene and polyethylene	Black	305 x 127 x 29.5



Staining jar, heavy glass, Coplin

Holds 10 microscope slides 76mm x 26mm back to back. Heavy glass with drop-on lid.

Catalogue No	Legacy No	
12858735	-	



Staining jar, polypropylene, Coplin



Scienceware®
BEL-ART PRODUCTS

This polypropylene Coplin staining jar improves on the classic tapered square design with a larger opening for easier placement and removal of slides.

The polypropylene screw cap is 53mm in diameter, has a flat top, easy grip larger knurled edge, deep thread and a seal ring to prevent leaks. Five internal slots hold up to 10 standard 76mm x 26mm slides back to back. The base diameter is 59.7mm, height is 114mm. Steam autoclavable at 121°C. Remove cap when autoclaving.

Catalogue No	Legacy No	Alt. No
10082844	-	F44208-0000



Staining set, microscope slide



Kartell

The trough is supplied with two lids; one lid for when the rack is not in use to prevent evaporation of staining liquid and the second with a 'cut out' to accommodate the rack handle. Made from highly resistant PMP. The polypropylene rack which holds 20 slides, 76mm x 26mm. Can be used with other troughs.

Trough dimensions (w x d x h), 81mm x 101mm x 86.5mm.

Rack dimensions (w x d x h), 69.5mm x 86mm x 21mm.

Not to be used with xylene.

Catalogue No	Legacy No	Alt. No
12382168	-	353X



All products are for research use only, unless otherwise indicated.