

Petri dishes, polystyrene, Thermo Scientific Sterilin



Single use, with an excellent uniform finish.

The dishes fit popular makes of automatic agar filling machines and the lids are recessed for stable stacking. CE marked in accordance with the In Vitro Diagnostic Device Directive 98/79/EC. 'Biologically clean'. Sterile by aseptic production.

Standard dishes



Catalogue No	Legacy No	Alt. No	Type	Dish diameter, mm	Height, mm	Pack qty
10655821	-	124	Deep, single vent	50	20.3	500
11359273	-	122	Single vent	50	14.5	700
11389273	-	101R20	Single vent	90	16.0	500
11309283	-	101VR20	Triple vent	90	16.2	500
11339283	-	501V	Triple vent	140	21.1	80
11379273	-	504	Contact plate	55	10.4	300
11329283	-	121V	Shallow, single vent	30	11.0	800
11369273	-	123	Single vent	60	15.1	540

Coloured dishes

Catalogue No	Legacy No	Alt. No	Type	Dish diameter, mm	Height, mm	Pack qty
11399273	-	101VBUE	Triple vent, blue	90	16.2	500
11319283	-	101VRED	Triple vent, red	90	16.2	500

Petri dishes, 100mm square (†lid turns 90° to allow venting or non-venting of the dish)

Catalogue No	Legacy No	Alt. No	Type	Dimensions, mm	Height, mm	Pack qty
11349273	-	109	Non-vented	100 x 100	21.1	120
11339273	-	103	25 compartments.†	100 x 100	20.8	120



All products are for research use only, unless otherwise indicated.

Dishes, polystyrene, with lids, Thermo Scientific Nunc® Nuncion™ Δ



- Cell culture treated
- Optically clear and very flat for use in microscopy
- 40mm and 58mm dishes available with grids for cloning or determination of plating efficiency
- Non-toxic
- Certified surface treatment for optimal cell attachment and growth

10508921, 10049030 and 10098720 are supplied in resealable bags.
10059030 are square dishes.

Catalogue No	Legacy No	Alt. No	Grids	Culture area, cm ²	Airvent	Dimensions, external, mm	Volume, working, mL	Inner pack qty	Pack qty
10654161	-	150318	-	8.8	No	40 x 12	3	10	500
10738581	-	153066	-	8.8	Yes	40 x 12	3	10	500
10191341	-	174926	2 x 2	8.8	Yes	40 x 12	3	20	500
10120101	-	150326	-	21.5	No	58 x 15	5	10	400
10111351	-	150288	-	21.5	Yes	58 x 15	5	10	400
10508921	-	150350	-	56.7	Yes	92 x 17	12.5	10	150
10049030	-	172958	-	56.7	Yes	92 x 21	12.5	10	480
10098720	-	168381	-	145	Yes	144 x 21	35	10	80
10059030	-	166508	-	500	Yes	245 x 245 x 25	135	4	16

