

Pipettes, Pasteur, polyethylene, transfer, Thermo Scientific



Micro

Catalogue No	Legacy No	Alt. No	Description	Length, mm	Bulb draw, mL	Drops, mL	Drops size	Graduations, mL	Pack qty
11337863	-	LW4237	Extended fine tip, sterile	153	3.0	50	20	-	500
11621839	-	LW4020	Standard	150	3.5	25	40	-	500
11377863	-	LW4252	Large bulb	150	7.7	22	45	-	250
11397863	-	LW4243	Fine tip ultra	50	1.0	45	22	-	500

Small

Catalogue No	Legacy No	Alt. No	Description	Length, mm	Bulb draw, mL	Drops, mL	Drops size	Graduations, mL	Pack qty
12334498	-	LW4273	3 inch	87	1.0	25	40	-	500
12344498	-	LW4295	Paddle	124	0.2	20	50	-	500
11307873	-	LW4293	Fine tip paddle	130	0.2	33	30	-	500

Long form

Catalogue No	Legacy No	Alt. No	Description	Length, mm	Bulb draw, mL	Drops, mL	Drops size	Graduations, mL	Pack qty
11641839	-	LW4694	Long reach	155	3.0	-	-	-	500
10111360	-	LW4050	9 inch	230	3.5	25	40	-	400
11317873	-	LW4051	9 inch sterile 1's	230	3.5	25	40	-	400
11327873	-	LW4052	9 inch sterile 10's	230	3.5	25	40	-	400
10705502	-	LW4053	9 inch sterile 20's	230	3.5	25	40	-	400
11651839	-	LW4263	Jumbo	300	10.0	22	45	-	100
11337873	-	LW4264	Jumbo sterile 1's	300	10.0	22	45	-	100
11347873	-	LW4266	Jumbo sterile 20's (pack of 100)	300	10.0	22	45	-	100

Standard

Catalogue No	Legacy No	Alt. No	Description	Length, mm	Bulb draw, mL	Drops, mL	Drops size	Graduations, mL	Pack qty
11357873	-	LW4336	Medi	155	2.4	25	40	-	500
11661839	-	LW4000	Standard	150	3.5	22	45	-	500
11367873	-	LW4003	Standard sterile 1's	150	3.5	22	45	-	500
11387873	-	LW4002	Standard sterile 20's	150	3.5	22	45	-	500
12394508	-	LW4204	Maxi	150	6.0	22	45	-	400
11397873	-	LW4205	Maxi, sterile 1's	150	6.0	22	45	-	400
11307883	-	LW4206	Maxi, sterile 20's	150	6.0	22	45	-	400
11317883	-	LW4141	Mayday	156	6.0	25	40	-	400
10793982	-	LW4691	Cell saver	137	4.6	N/A	62	-	250

Pipettes, Pasteur, polyethylene, Thermo Scientific Sterilin®



One piece polyethylene moulded pipette. Maximum bulb capacity approx. 3.5mL. Length 155mm.

Catalogue No	Legacy No	Alt. No	Description	Capacity, mL	Inner pack qty	Pack qty
12999477	-	202C	Fine bore, non-sterile	1	-	500
11329203	-	201C	Non-sterile, graduated	1	-	500
11339203	-	PP88SA	Sterile, graduated	1	10	500
11349203	-	PP88SB	Sterile, graduated	1	1	500
12909487	-	200C	Non-sterile, graduated	3	-	500
11369203	-	PP89SA	Sterile, graduated	3	10	500
11379203	-	PP89SB	Sterile, graduated	3	1	500

