

Pipettor tips, speciality, universal fit, Fisherbrand® SureOne™



Fisherbrand®

Fisherbrand® SureOne™ speciality pipettor tips include gel loading, genomic and extended length tips, in bulk, racked and sterile racked packaging.

SureOne™ universal fit, speciality pipettor tips are compatible with most popular brands of pipettor, optimal fit with Fisherbrand® Elite™ pipettor.

Non-sterile products - certified to be free of RNase/DNase and DNA.

Sterile products - e-beam sterile products are certified to be free of RNase/DNase, pyrogen, bioburden, PCR* inhibitors and endotoxins.

200µL gel loading, outer diameter 0.58mm

Catalogue No	Legacy No	Volume, µL	Pack type	Colour	Inner pack qty	Pack qty
12962571	-	1 to 200	Bulk	Clear	-	1,000
12972571	-	1 to 200	Racked	Clear	96	960
12982571	-	1 to 200	Racked sterile	Clear	96	960

20µL gel loading, filter tip, outer diameter 0.58mm

Catalogue No	Legacy No	Volume, µL	Pack type	Colour	Inner pack qty	Pack qty
12992571	-	2 to 20	Filtered sterile	Clear	204	1,020

100µL gel loading, filter tip, outer diameter 0.58mm

Catalogue No	Legacy No	Volume, µL	Pack type	Colour	Inner pack qty	Pack qty
12902581	-	10 to 100	Filtered sterile	Clear	204	1,020

200µL extended length, 90mm long

Catalogue No	Legacy No	Volume, µL	Pack type	Colour	Inner pack qty	Pack qty
12912581	-	1 to 200	Bulk	Clear	-	1,000
12922581	-	1 to 200	Racked	Clear	204	1,632
12932581	-	1 to 200	Racked sterile	Clear	204	1,632

200µL extended length, filter tip, 90mm long

Catalogue No	Legacy No	Volume, µL	Pack type	Colour	Inner pack qty	Pack qty
12942581	-	20 to 200	Filtered sterile	Clear	204	1,632

200µL large orifice tip, inner diameter 1.5mm

Catalogue No	Legacy No	Volume, µL	Pack type	Colour	Inner pack qty	Pack qty
12952581	-	1 to 200	Bulk	Clear	-	1,000
12962581	-	1 to 200	Racked	Clear	96	960
12972581	-	1 to 200	Racked sterile	Clear	96	960

200µL large orifice filter tip, inner diameter 1.5mm

Catalogue No	Legacy No	Volume, µL	Pack type	Colour	Inner pack qty	Pack qty
12982581	-	20 to 200	Filtered sterile	Clear	96	960

1,000µL large orifice tip, graduated at 100µL, 200µL, 500µL and 1,000µL, inner diameter 1.5mm

Catalogue No	Legacy No	Volume, µL	Pack type	Colour	Inner pack qty	Pack qty
12992591	-	100 to 1,000	Bulk	Clear	-	1,000
12902591	-	100 to 1,000	Racked	Clear	96	960
12912591	-	100 to 1,000	Racked sterile	Clear	96	960

1,000µL large orifice filter tip, graduated at 100µL, 200µL, 500µL and 1,000µL, inner diameter 1.5mm

Catalogue No	Legacy No	Volume, µL	Pack type	Colour	Inner pack qty	Pack qty
12922591	-	100 to 1,000	Filtered sterile	Clear	96	960



*Polymerase Chain Reaction (PCR) is a process covered by patents owned by Hoffmann-La Roche