



Pipettor tips, Thermo Scientific Finn timer®



Accuracy and precision depend on the quality of the tips used for pipetting. Designed to optimise the performance of Thermo Scientific Finn timer®, advanced moulding techniques and strict quality control ensure a perfect fit. Specialist tips, including the Finn timer® 250 and 1000 range, have a bevelled orifice which reduces shear force and thereby gives added security when pipetting fragile items.

Finn timer® 10/20/50 are autoclavable.

Finn timer® 10

Catalogue No	Legacy No	Alt. No	Volume range, µL	Pack type	Inner pack qty	Pack qty
10666491	-	9400310	0.2 to 10	Bag	-	1,000
10676491	-	9400300	0.2 to 10	Rack	96	960
10165280	-	9400303	0.2 to 10	Sterile rack	96	960
12937047	-	9400326	0.2 to 10	Rack refill starter kit	192	2,688
10633323	-	9400327	0.2 to 10	Rack refill kit	192	3,840

Finn timer® Flex™ 10

Catalogue No	Legacy No	Alt. No	Volume range, µL	Pack type	Inner pack qty	Pack qty
11863420	-	94060100	0.1 to 10	Rack	96	960
11873420	-	94060103	0.1 to 10	Sterile rack	96	960
11873430	-	94060116	0.1 to 10	Refill starter kit	96	2,016
11883430	-	94060117	0.1 to 10	Refill kit	96	1,920
11883420	-	94060120	0.1 to 10	Bag	-	1,000

Finn timer® 20

Catalogue No	Legacy No	Alt. No	Volume range, µL	Pack type	Inner pack qty	Pack qty
10061292	-	9400610	0.2 to 20	Rack	384	3,840
10132582	-	9400613	0.2 to 20	Sterile rack	384	3,840
10152532	-	9400620	0.2 to 20	Bag	-	1,000

Finn timer® 50

Catalogue No	Legacy No	Alt. No	Volume range, µL	Pack type	Inner pack qty	Pack qty
10423322	-	9400370	0.2 to 50	Rack	384	3,840
10799464	-	9400373	0.2 to 50	Sterile rack	384	3,840
10512861	-	9400360	0.2 to 50	Bag	-	1,000

Finn timer® wide orifice, 10µL to 250µL

Catalogue No	Legacy No	Alt. No	Volume range, µL	Pack type	Inner pack qty	Pack qty
10389220	-	9405020	10 to 250	Bag	-	1,000
10165710	-	9405030	10 to 250	Bulk	-	22,000
10247330	-	9405120	10 to 250	Rack	96	960
10349320	-	9405123	10 to 250	Sterile rack	96	960

Finn timer® 200, extended, 78mm length

Catalogue No	Legacy No	Alt. No	Volume range, µL	Pack type	Inner pack qty	Pack qty
10348890	-	9400100	5 to 200	Box	-	400
10503241	-	9400130	5 to 200	Rack	96	960
10607841	-	9400133	5 to 200	Sterile rack	96	960

Finn timer® Flex™ 200

Catalogue No	Legacy No	Alt. No	Volume range, µL	Pack type	Inner pack qty	Pack qty
10411844	-	94060310	1 to 200	Rack	96	960
10335484	-	94060313	1 to 200	Sterile rack	96	960
11893430	-	94060316	1 to 200	Refill starter kit	96	2,016
11803440	-	94060317	1 to 200	Refill kit	96	1,920
11813430	-	94060320	1 to 200	Bag	-	1,000