



Pipettors, dispensing, Thermo Scientific FinnpiPETTE® Stepper, Multistepper and Stepper tips



- Makes up to 45 deliveries before refilling is necessary
- Positive displacement design ensures accurate results, even for viscous, volatile or other problematic liquids
- Lightweight, ergonomic handle and wide finger rest minimise operator hand strain
- Operation is simple, just attach the appropriate size tip, fill and prime the tip and dial the selected volume
- Volume is easy to set, with the aid of the convenient volume setting chart located on the handle

Ordering alerts: Universal stepper handle is also available with an eight channel module.

Required accessories: Uses the Thermo Scientific FinnpiPETTE® Stepper tip, a plunger style tip available in sterile or non-sterile form. This tip comes in seven sizes from 0.5mL to 50mL.

Catalogue No	Legacy No	Alt. No	Description	Volume range, µL
10358980	-	4540000	FinnpiPETTE® Stepper, single channel, handle	10 to 5,000
10155280	-	4540500	FinnpiPETTE® Multistepper, 8 channel, handle and module	50 to 250
10591151	-	2206590	FinnpiPETTE® Multistepper, 8 channel, module only	50 to 250

Stepper tips Non-sterile

Catalogue No	Legacy No	Alt. No	Tip size	Dispensing volume, µL	Pack qty
10626111	-	9404170	0.5mL	10, 20, 30, 40, 50	100
10185900	-	9404180	1.25mL	25, 50, 75, 100, 125	100
10542481	-	9404190	2.5mL	50, 100, 150, 200, 250	100
10762111	-	9404200	5.0mL	100, 200, 300, 400, 500	50
10451311	-	9404210	12.5mL	250, 500, 750, 1,000, 1,250	50
10401891	-	9404220	25mL	500, 1,000, 1,500, 2,000, 2,500	20
10083630	-	9404230	50mL	1,000, 2,000, 3,000, 4,000, 5,000	10

Sterile, individually packed

Catalogue No	Legacy No	Alt. No	Tip size	Dispensing volume, µL	Pack qty
10411891	-	9404173	0.5mL	10, 20, 30, 40, 50	50
10175280	-	9404183	1.25mL	25, 50, 75, 100, 125	50
10237620	-	9404193	2.5mL	50, 100, 150, 200, 250	50
10637071	-	9404203	5.0mL	100, 200, 300, 400, 500	25
10740791	-	9404213	12.5mL	250, 500, 750, 1,000, 1,250	25
10023660	-	9404223	25mL	500, 1,000, 1,500, 2,000, 2,500	10
10750981	-	9404233	50mL	1,000, 2,000, 3,000, 4,000, 5,000	10

Multistepper tips

Catalogue No	Legacy No	Alt. No	Description	Inner pack qty	Pack qty
10033680	-	9401300	Box	-	400
10145520	-	9401330	Rack	96	960
10155520	-	9401333	Sterile rack	96	960
10514561	-	94052450	Filter sterile rack	96	960

Accessories

Catalogue No	Legacy No	Alt. No	Description	Pack qty
10665911	-	9420 300	Adapter	10
10606111	-	9420310	Adapter sterile	5
10481121	-	9510027	Reagent reservoir, 'V' shaped, autoclavable	5
10207340	-	9510037	8 channel reservoir, autoclavable for 4, 8 or 12 channel pipettors	10
10491121	-	9510047	Single channel reservoir, autoclavable	10