

Spatulas, vibrating, Quaver



Make accurate additions of precious materials as needed. Spinning the toothed thumb wheel which vibrates the spatula, gives control of how much or how little is used. Available with a plain stainless steel blade or with a brown Fluo-Kem® PTFE coating containing FEP. The coating complies with FDA regulations governing direct food contact applications and reduces sticking losses. 180mm long.



Catalogue No	Legacy No	Alt. No	Description
12388109	-	H36714-0000	Stainless steel blade
11307184	-	H36715-0000	Brown FTPE FEP coated stainless steel blade

Spatulas, vibrating, Quaverette



Offers the same features as the Quaver vibrating spatula, but permits utilisation of smaller quantities and containers. The narrow v-shaped blade and 210mm length manoeuvres easily into small mouthed containers and tight corners. Available with brown Fluo-Kem® coating containing FEP which complies with FDA regulations governing direct food contact applications.



Catalogue No	Legacy No	Alt. No	Description
11362204	-	H36718-0000	Stainless steel blade
11317184	-	H36719-0000	Brown Fluo-Kern® PTFE containing FEP coated stainless steel blade

Spatula weighing set, micro



Set of 6 micro size double-ended tools and a pair of tweezers for accurate handling of minute quantities.

Each tool has two different end shapes, doubling the selection of implements to better match the users' needs. Tool ends include two small spoons for transferring and collecting both dry and liquid materials. There are also eight flat spatulas of various widths, lengths and end shapes, that are ideal for mixing, scraping and spreading. The tweezers have a serrated inside surface for secure sample handling and grooved handles for better gripping. Approximate length of instruments 140mm.

All tools are manufactured from polished stainless steel and can be steam autoclaved at 121°C. When not in use the tools can be stored inside the carrying case provided.



Catalogue No	Legacy No	Alt. No	
11352204	-	F36705-0000	

Spatula, wooden applicator

Catalogue No	Legacy No	Pack qty
12682560	-	864



For scoops, please see page 1366

