

### Stirrer hotplates, C-Mag



### Technical Support

For further information on any of the featured products contact your local office



- Ceramic glass plate, offering excellent chemical resistance
- Fixed safety circuit of 550°C
- Hot surface warning indicator to prevent burns
- Stirring quality 5L, 10L and 15L (H<sub>2</sub>O)
- Digital display
- Digital error display
- Elevated control panel for protection against leaking liquids

#### Technical Specification - Specific

| Model No                              | HS 4            | HS 7            | HS 10           |
|---------------------------------------|-----------------|-----------------|-----------------|
| Plate dimensions, mm                  | 100 x 100       | 180 x 180       | 260 x 260       |
| Stirring quantity H <sub>2</sub> O, L | 5               | 10              | 15              |
| Dimensions [w x d x h], mm            | 150 x 260 x 105 | 220 x 330 x 105 | 300 x 415 x 105 |
| Output                                | 250W (heat)     | 1,000W (heat)   | 1,500W (heat)   |
| Mass, kg                              | 3               | 5               | 6               |

#### Technical Specification - General

|                       |               |
|-----------------------|---------------|
| Temperature range, °C | 50 to 500     |
| Speed range, rpm      | 100 to 1,500  |
| Display               | Digital       |
| Electrical supply     | 230V, 50/60Hz |

| Catalogue No | Legacy No | Alt. No | Description                                                                                               |
|--------------|-----------|---------|-----------------------------------------------------------------------------------------------------------|
| 10367894     | -         | 3581002 | C-MAG HS 4 hotplate stirrer                                                                               |
| 10625435     | -         | 3581202 | C-MAG HS 7 hotplate stirrer                                                                               |
| 10347904     | -         | 3581402 | C-MAG HS 10 hotplate stirrer                                                                              |
| 12917246     | -         | 9015902 | C-MAG HS 7 hotplate package c/w hotplate stirrer, ETS-D5 contact thermometer and support and holding rods |

#### Accessories

| Catalogue No | Legacy No | Alt. No | Description                                 |
|--------------|-----------|---------|---------------------------------------------|
| 10555243     | -         | ETS-D5  | Contact thermometer (for HS7 and HS10 only) |
| 10120881     | -         | 1545100 | H16V support rod                            |
| 10043852     | -         | 3547700 | H38 Holding rod                             |
| 10058550     | -         | 2437700 | H44 Boss head                               |
| 10321352     | -         | 1293100 | Stir bar remover                            |