

Whether you require glassware or digital thermometers, we supply an extensive range featuring the industry's leading brands. Here you will find a wide selection of interchangeable probes including infrared, rigid, flexible, air and penetration models to ensure precision in the lab whatever the application. Printing/logging, alarms and waterproof models are also available, providing choice and accuracy at a great price.



Thermometers, white backed, red spirit filled, total immersion, Fisherbrand®

Fisherbrand®



General purpose/stirring thermometers with amber graduations, reinforced bulb for stirring. Non-roll fitting.

Catalogue No	Legacy No	Range, °C	Subdivisions, °C	Length, mm
11577273	-	-10 to 50	0.5	305
11587273	-	-20 to 110	1	305
11597273	-	-20 to 150	1	305
11507283	-	-10 to 50	1	155
11517283	-	-10 to 110	1	155

Thermometers, white backed, red spirit filled, 76mm immersion, Fisherbrand®

Fisherbrand®



General purpose/stirring thermometer with amber graduations, reinforced bulb for stirring. Non-roll fitting.

Catalogue No	Legacy No	Range, °C	Subdivisions, °C	Length, mm
11551993	-	-10 to 50	0.5	305
11561993	-	-10 to 110	1	305
11571993	-	-10 to 150	1	305