



Vials, clear tubular glass, 7mL Bijou and 28mL, universal, with fitted polypropylene screw caps, Fisherbrand®



Type III clear soda glass vials with screw neck and fitted PP cap

Catalogue No	Legacy No	Height, mm	Diameter, mm	Capacity, mL	Tray quantity	Pack qty
14803562	-	43	20	7	1 x tray 245	245
14813562	-	43	20	7	7 x trays 245	1,715
14823562	-	82	27	28	1 x tray 125	125
14833562	-	82	27	28	4 x trays 125	500



Vials, clear moulded glass, pathology media, with fitted screw caps, Fisherbrand®



Type I clear neutral glass vials with screw neck and fitted PP cap

Catalogue No	Legacy No	Height, mm	Diameter, mm	Capacity, mL	Pack qty
14843562	-	50	22	7	288
14863562	-	83	28	28	144

Type I clear neutral glass vials with screw neck and fitted aluminium rubber lined caps

Catalogue No	Legacy No	Height, mm	Diameter, mm	Capacity, mL	Pack qty
14873562	-	50	22	7	20
14853562	-	50	22	7	288
14883562	-	65	28	14	20
14893562	-	65	28	14	288
14803572	-	83	28	28	20
14823572	-	83	28	28	144



Vials, clear neutral glass, Type I, with dropper assembly, Fisherbrand®



- Clear Type I neutral glass vials with screw neck
- Supplied with dropper assembly comprising clear glass dropper tube fitted to a polypropylene cap with natural rubber bulb
- Droppers are assembled but not fitted to vials (supplied in a separate bag)
- Dropper assemblies for 3mL and 4.5mL vials have white caps and bulbs, 5mL and 10mL vials have black caps and bulbs

Catalogue No	Legacy No	Height, mm	Diameter, mm	Capacity	Neck diameter, mm	Pack qty
14813572	-	36	16	3mL vial, 40µL dropper	14.5	50
14833572	-	36	16	3mL vial, 40µL dropper	14.5	374
14843572	-	47	16	4.5mL vial, 40µL dropper	14.5	50
14853572	-	47	16	4.5mL vial, 40µL dropper	14.5	374
14863572	-	43	20	5mL vial, 50µL dropper	18	40
14873572	-	43	20	5mL vial, 50µL dropper	18	245
14883572	-	58	20	10mL vial, 50µL dropper	18	40
14893572	-	58	20	10mL vial, 50µL dropper	18	245

For Fisherbrand® colour coder vial caps, please see page 1602.

