



Contamination Control

Wipes



Wipes, cleanroom, ISO Class 6-7

04



NEW

TechniCloth offers the lowest level of releasable ions and non-volatile contaminants, contains no chemical binders and tolerates most standard solvents. Its blend of 55% cellulose and 45% polyester provides the optimum balance of sorption, durability, and cleanliness.

Catalogue No	Alt. No	Dimensions, mm	Pack qty
CMH-820-030V	TX609	230 x 230	300
CMH-820-050P	TX612	310 x 310	150



Wipes, cleanroom, Class 100, ISO Class 5, Crew 2®

04



Single ply non-woven fabric, bonded to eliminate and trap surface particles. Low lint and extractors. Double packaged in anti-static polybags.

Catalogue No	Alt. No		Pack qty
CMH-630-011L	7605		500

Stop, think...

have you got the right hand protection?

Comprehensive range of PPE and Safety Products, we make your safety our business

