



Beakers, Pyrex® glass, squat form, ISO3819

- With spout and graduations except 10,000mL size

Catalogue No	Cap., mL	Height, mm	O.D., mm	Pack qty
BNB-320-010N	5	30	22	10
BNB-320-030H	10	35	26	10
BNB-320-050B	25	50	34	10
BNB-320-070S	50	60	42	10
BNB-320-090M	100	70	50	10
BNB-320-110J	150	80	60	10
BNB-320-130D	250	95	70	10
BNB-320-150U	400	110	80	10
BNB-320-190Y	600	125	90	10
BNB-320-215S	800	135	100	10
BNB-320-235M	1,000	143	105	10
BNB-320-275A	2,000	185	130	10
BNB-320-295R	3,000	210	150	4
BNB-320-310B	5,000	270	170	each
BNB-320-330S	10,000	350	217	each

Beakers, Pyrex® glass, heavy duty, squat form

- With spout and double graduated metric scale

Catalogue No	Cap., mL	Height, mm	O.D., mm	Pack qty
BNB-400-020Q	150	89	57	12
BNB-400-031L	250	90	68	12
BNB-400-051F	400	110	77	12
BNB-400-071W	600	124	90	6
BNB-400-111N	1,000	156	108	6
BNB-400-131H	2,000	179	131	4
BNB-400-140G	3,000	230	170	each
BNB-400-150D	4,000	250	160	each
BNB-400-160A	5,000	300	180	each



Beakers, ultra clear polypropylene, **Fisherbrand**® squat form, BS5404 Pt1, ISO7056



- Large pouring spout and moulded graduations
- Not suitable for stirrer hotplates

Catalogue No	Cap., mL	Pack qty
FB55910	25	10
FB55912	50	10
FB55914	100	10
FB55916	250	10
FB55918	500	10
FB55920	1,000	5

Beakers, PTFE, squat form



- Isostatically moulded from pure PTFE
- Non-contaminating and completely inert to highly corrosive reagents
- Non-wettable, lightweight and heat resistant
- Spout for easy pouring
- Can be heated on a hotplate up to 280°C when proper care is exercised

Catalogue No	Cap., mL
FB59989	5
FB59990	10
FB59991	25
FB59992	50
FB59993	100
FB59994	250
FB59995	500
FB59996	1,000
FB59997	2,000
FB59998	5,000

