

Petri dishes



Petri dishes, polystyrene

Fisherbrand®

- Manufactured from very clear crystal polystyrene, extremely flat with excellent optical qualities and non-distorting below 75°C
- The melting temperature of the raw material is approximately 250°C (This is considered to be sufficient to destroy bacteria)
- The dishes are packed under very strict aseptic conditions

Catalogue No	Dia., mm	Type	Sterile	Pack qty
FB51480	55	Triple vent	Yes	1,620
FB51483	140	Triple vent	Yes	176
FB51502	55	Triple vent	No	1,620
FB51504	90	Single vent	No	600
FB51506	90	Triple vent	No	600
FB51508	90	Single vent	Yes	600
FB51510	90	Triple vent	Yes	600



Petri dishes, polystyrene

- Single use, with an excellent uniform finish
- Fit popular makes of automatic agar filling machines
- Recessed, for stable stacking
- Produced and packed under aseptic conditions

Catalogue No	Type	Dia., mm	Pack qty
PDS-140-020X	Single vent	50	700
PDS-140-045V	Single vent	90	500
PDS-140-050F	Triple vent	90	500
PDS-140-055S	2 compartment	90	500
PDS-140-060C	3 compartment	90	500
YXA-530-060N	3 compartment	90	20
PDS-140-080T	Triple vent	140	80
PDS-140-085J	Contact plate	55	300
PDS-140-120K	Square	100 x 100	120
PDS-140-130H	Square Repli., 25 compartments	100 x 100	120
PDS-140-140E	Triple vent, amber	90	500
PDS-140-150F	Triple vent, blue	90	500
PDS-140-160V	Triple vent, red	90	500
PDS-140-015H	Deep, single vent	50	500
PDS-140-170S	Shallow, single vent	30	800
PDS-140-180P	Single vent	60	540



Petri dishes, Poxygrid® rack



- Epoxy coated steel construction which holds 100mm diameter petri dishes
- Ideal for lab transport and the incubation of dishes
- Holds 6 x 100mm dishes
- Dimensions (w x d x h), mm: 130 x 85 x 120

Catalogue No
PDS-270-010L