

page
375

Evaporating basins, porcelain,
round bottom with spout



- Shallow form
- Glazed inside and outside

Catalogue No	Cap., mL	Height, mm	Dia., mm
BHW-380-030F	50	30	75
BHW-380-050W	75	32	80
BHW-380-070Q	95	36	85
BHW-380-090K	115	40	100

beakers

page
383

Beakers, Pyrex® glass, squat
form, ISO3819

PYREX®



- With spout and graduations except 10,000mL size

Catalogue No	Cap., mL	Height, mm	O.D., mm	Pack qty
BNB-320-070S	50	60	42	10
BNB-320-090M	100	70	50	10
BNB-320-110J	150	80	60	10
BNB-320-130D	250	95	70	10
BNB-320-150U	400	110	80	10
BNB-320-190Y	600	125	90	10

basins

page
382

Beakers, borosilicate glass, squat form

Fisherbrand®



- With spout and graduations

Catalogue No	Cap., mL	Height, mm	Dia., (external), mm	Pack qty
FB33109	50	60	42	10
FB33110	100	70	50	10
FB33111	150	80	60	10
FB33112	250	95	70	10
FB33113	400	110	80	10
FB33116	1,000	143	105	each
FB33117	2,000	185	130	each

page
375

Beakers, polypropylene,
squat form, BS5404 : Part 1, ISO7056

Fisherbrand®



- Translucent with easy to read moulded graduations
- Large pouring spout
- Not suitable for stirrer hot plates

Catalogue No	Cap., mL	Pack qty
FB55910	25	10
FB55912	50	10
FB55914	100	10
FB55916	250	10
FB55918	500	10
FB55920	1,000	5