

## Touch-Counter™

### Cuts counting time in half

Touch-Counter™ is a combination marker and digital counter designed for the lab and plant. Each time an item is touched, it is automatically marked by the felt-tip pen and counted by the built-in digital counter. Use it to count sample bottles, cell growths, and to receive/inventory/ship supplies.



### One-touch operation saves time

Operation is simple: Touch an item and it is marked and counted simultaneously. Impossible to make an error, even when interrupted. Marks hold your place and counter holds your total. Durable aluminium unit has an audible click, zero reset knob, and a 9999-digit counter. Touch-Counter™ is 178 x 51 x 38 mm in length and weight is 186 g. Spring-actuated ink reservoir valve seals felt tip when not in use and prevents ink dry-out.

| Description   | Cat. No.       | Price   |
|---|----------------|---------|
| Fisher Scientific Touch-Counter™—<br>Supplied with black disposable pen.                      | <b>FB70281</b> | £105.00 |
| Disposable 5-Pen Pack™—<br>Package contains one each: black, red, blue, green,<br>and yellow. | <b>FB70282</b> | £50.00  |

## Counter-Pen™

Counter-Pen™ is a combination marker and digital counter designed for the lab and plant. Each time a black mark is made with the felt-tipped pen, it automatically sounds a beep, and LCD displays the count. Use it to count cell growths on petri dishes, inventory parts, and tally samples received.

Touch, sound, and count occur simultaneously. Impossible to make an error, even when interrupted—marks hold your place and counter holds your total. Durable ABS plastic unit provides audio/visual counts to 99999 with reset key, on/off switch, and pen cap. Size is 152 x 17 mm and weight is 31 g.

| Description                    | Cat. No.       | Price  |
|--------------------------------|----------------|--------|
| Fisher Scientific Counter-Pen™ | <b>FB70375</b> | £20.00 |



## Digital Key-Chain Counter/ Wrist Strap

Key-Chain/Wrist-Strap/Carabiner Digital Counter is designed for ease-of-use and fail-safe counting. Features include display count to 99,999, two year constant use battery, water-resistant O-ring-sealed with shock-resistant ABS plastic case, key click for count confirmation, and recessed no accident reset.

Eliminates gear and miscount problems associated with mechanical counters. Diameter is 47.6 x 12.7 mm weight is less than 29 g.



## Digital Counter

Lightweight counter reads to 99,999. Water-resistant, O-ring-sealed and shock-resistant ABS plastic case is ideal for lab and plant use. Tactile click keys provide positive action.

### Completely electronic

Easy-to-read LCD and solid-state construction eliminate gear and miscount problems associated with mechanical counters. Fail-safe design eliminates accidentally resetting. Size is 67 x 42 x 21 mm and weight is 31 g. Battery powers unit for over 2 years of continuous operation.



| Description                       | Cat. No.       | Price  |
|-----------------------------------|----------------|--------|
| Fisher Scientific Digital Counter | <b>FB70280</b> | £24.50 |

| Description   | Cat. No.       | Price  |
|---|----------------|--------|
| Fisher Scientific Digital Key-Chain Counter/Wrist Strap | <b>FB61154</b> | £11.50 |