




Petri dishes, polystyrene



56

Single use, with an excellent uniform finish.

The dishes fit popular makes of automatic agar filling machines and the lids are recessed for stable stacking. CE marked in accordance with the In Vitro Diagnostic Device Directive 98/79/EC. 'Biologically clean'. Sterile by aseptic production.

Standard dishes

Catalogue No	Alt. No	Type	Dish diameter, mm	Height, mm	Pack qty
PDS-140-015H	124	Single vent, deep	50	20.3	500
PDS-140-020X	122	Single vent	50	14.5	700
PDS-140-045V	101R20	Single vent	90	16.0	500
PDS-140-050F	101VR20	Triple vent	90	16.2	500
PDS-140-080T	501V	Triple vent	140	21.1	80
PDS-140-085J	504	Contact plate, triple vent	55	10.4	300
PDS-140-170S	121V	Triple vent	30	11.0	800
PDS-140-180P	123	Single vent	60	15.1	540

Coloured dishes

Catalogue No	Alt. No	Type	Dish diameter, mm	Height, mm	Pack qty
PDS-140-140E	101VAMB	Triple vent, amber	90	16.2	500
PDS-140-150F	101VBLUE	Triple vent, blue	90	16.2	500
PDS-140-160V	101VRED	Triple vent, red	90	16.2	500

Petri dishes, 100mm square (lid turns 90° to allow venting or non-venting of the dish)

Catalogue No	Alt. No	Type	Dimensions, mm	Height, mm	Pack qty
PDS-140-120K	109	Square, non-vented	100 x 100	21.1	120
PDS-140-130H	103	Square, 25 compartments	100 x 100	20.8	120

Petri dishes, polystyrene



76

For culturing of fungi, bacteria and other microorganisms.

- Perform well in automatic dispensers due to complete flatness and uniform height
- Deep petri dishes allow for longer culture periods
- Compatible with automated systems

Catalogue No	Alt. No	Dimensions, internal [dia x d], mm	Dimensions, external, mm	Sterile or non-sterile	Culture area, cm ²	Total volume, mL	Inner pack qty	Pack qty
DIS-985-030B	240401	136 x 18	145 x 21 (includes lid)	Non-sterile	145	250	10	180
TKT-310-010C	254925	86 x 12	100 x 15 (includes lid)	Non-sterile	58	68	20	320
TKT-310-020W	263991	86 x 12	100 x 15 (includes lid)	Sterile	58	68	20	320
TKT-310-450U	249964	136 x 18	145 x 20 (includes lid)	Sterile	145	250	10	80

All products are for research use only unless otherwise indicated.