

Cryogenic vial inserts, colour coded



55

Provide quick visual identification of specific vials or type of sample.

Polystyrene discs which fit into the top of the cryogenic vial closure.

Cut-out tops to facilitate handling of vials with forceps and flat surfaces to permit writing on top of vials. Supplied in packs of 100.

Catalogue No	Alt. No	Colour
CRY-100-510N	5045-0000	White
CRY-100-512J	5045-0002	Yellow
CRY-100-514F	5045-0003	Blue
CRY-100-516B	5045-0004	Green
CRY-100-518U	5045-0005	Red



Cryogenic, vials, CryoELITE™




04

For storing, tracking and transporting samples.

- Conical well design allows precise sample retrieval for both automated and manual systems
- Free standing vials have Loctagon™ vial skirt which provides stability in the freestanding position
- Free standing vials lock into the CryoELITE benchmate rack to provide easy one handed manipulation
- Accommodate 2D data matrix bar code bottom insert to provide a unique identifier for traceability
- Certified RNase/DNase and Endotoxin Free
- Sterile
- Non-pyrogenic

Internal thread, freestanding

Catalogue No	Alt. No	Dimensions [dia. x h], mm	Capacity, mL	Cap colour	Writing area [Y/ N]	Sterile, Yes/No	Pack qty
CRY-161-005B	W985915	12 x 40	1.2	Natural	Y	Yes	500
CRY-161-010Y	W985902	12 x 50	2	Natural	Y	No	1000
CRY-161-015V	W985903	12 x 50	2	Natural	N	No	1000
CRY-161-020F	W985922	12 x 50	2	Natural	Y	Yes	500
CRY-161-025S	W985916	12 x 50	2	White	Y	Yes	500
CRY-161-030C	W985917	12 x 50	2	Red	Y	Yes	500
CRY-161-035P	W985918	12 x 50	2	Pink	Y	Yes	500
CRY-161-040W	W985919	12 x 50	2	Yellow	Y	Yes	500
CRY-161-045M	W985920	12 x 50	2	Green	Y	Yes	500
CRY-161-050T	W985921	12 x 50	2	Blue	Y	Yes	500
CRY-161-055J	W985923	12 x 63	3	Natural	Y	Yes	500
CRY-161-060Q	W985924	12 x 77	4	Natural	Y	Yes	500
CRY-161-065G	W985925	12 x 91	5	Natural	Y	Yes	500



Internal thread, round bottom. Cap colour: Natural, enhanced translucent.

Catalogue No	Alt. No	Dimensions [dia. x h], mm	Capacity, mL	Writing area [Y/ N]	Sterile, Yes/No	Pack qty
CRY-161-200D	W985910	12 x 39	1.2	Y	Yes	500
CRY-161-205Q	W985911	12 x 49	2	Y	Yes	500
CRY-161-210A	W985900	12 x 49	2	Y	No	1000
CRY-161-215N	W985901	12 x 49	2	N	No	1000
CRY-161-220U	W985912	12 x 63	3	Y	Yes	500
CRY-161-225K	W985913	12 x 76	4	Y	Yes	500
CRY-161-230R	W985914	12 x 90	5	Y	Yes	500



entry continued