

### Cryogenic vials, polypropylene, Cryule®



03

Externally threaded, these polypropylene vials and caps are freestanding on the bench and lock into the specially designed rack for one handed cap removal. The unique 50 position rack securely locks freestanding Cryule® vials in place and is made of durable autoclavable polypropylene and allows easy sample identification with an alpha-numeric index on top of each rack. Moulded corner posts provide stability when stacking (remove non-skid feet inserts). Non-skid feet offer additional stability for benchwork.

Note: Cryogenic vials to be used in gas phase of liquid nitrogen only.

Catalogue No	Alt. No	Cap colour	Diameter, mm	Height, mm	Pack qty
<b>CRY-160-030Y</b>	985734	White	12	46	500
<b>CRY-160-050C</b>	985735	Red	12	46	500
<b>CRY-160-070T</b>	985736	Pink	12	46	500
<b>CRY-160-090N</b>	985737	Yellow	12	46	500
<b>CRY-160-110K</b>	985738	Green	12	46	500
<b>CRY-160-130E</b>	985739	Blue	12	46	500

Accessory  
Rack

Catalogue No	Alt. No	Length, mm	Width, mm	Height, mm	Pack qty
<b>CRY-160-500U</b>	985810	190	190	20	5



### Cryogenic vials, Cryule®



03

Externally threaded. Uniquely designed vials are more compact, and unlike many other vials, safely store the total labelled volume. Supplied sterilised with caps attached, 50 per bag, 10 bags per case. White marking area. The 50-position rack is ideal for storing and transporting 1mL or 2mL Cryule® or even 12mm x 75mm glass culture tubes. Each rack is constructed of durable autoclavable polypropylene; moulded corner posts provide stability when stacking.

Catalogue No	Alt. No	Volume, mL	Diameter, mm	Height, mm	Pack qty
<b>CRY-170-030B</b>	985730	1	12	24	500
<b>CRY-170-050S</b>	985731	2	12	44	500

Accessory  
Cryule® polypropylene rack

Catalogue No	Alt. No	Width, mm	Depth, mm	Height, mm	Pack qty
<b>CRY-170-500N</b>	985800	190	100	20	5



### Cryogenic storage accessories, Cryo system



76

- Thermo Scientific Nunc CryoFlex™ is patented
- CryoTube™ racks are autoclavable
- Cans can be marked with CryoStore® tabs
- Nunc CryoFlex™ provides extra safety when freezing hazardous samples or storing in liquid nitrogen



Catalogue No	Alt. No	Description	Dimensions, mm	Material	Inner pack qty	Pack qty
<b>TES-950-525P</b>	347155	CryoTube® rack for 70 tubes	233 x 150	Polycarbonate	1	30
<b>CRY-960-501A</b>	376589	CryoTube® rack for 40 tubes	202 x 102	Polyphenyleneoxide	1	15
<b>CRY-960-503T</b>	375930	CryoColour® coders, assorted	-	Polystyrene	500	2,000
<b>CRY-960-505P</b>	343958	Cover tubing, Nunc CryoFlex™	500 x 14	Polyethylene	15	300
<b>CRY-960-507L</b>	378441	CryoStore® cane	-	Aluminium	1	50