


DNA Amplification

Plastics - PCR tubes & strips



*Polymerase Chain Reaction (PCR) is a process covered by patents owned by Hoffman-La Roche

Micro tubes and plates, thin wall, Purepak™



03

Efficient heat transfer

The thin wall construction ensures uniform contact with the thermal cycler block, producing optimal heat transfer. This reduces cycle times and helps maintain accurate reaction temperatures.

Reduced sample evaporation

Secure fitting caps that stay closed during thermal cycling but are still easy for the user to open. The secure seal also minimises sample evaporation.

Premium quality

Precision moulds, state of the art manufacturing equipment and premium-grade polypropylene. In addition, tubes are quality inspected and certified RNase, DNase, DNA and pyrogen free by independent labs. These extra measures result in consistently superior tubes that are unsurpassed in the industry. Non-sterile.

Purepak™

Each pack is divided into individually sealed bags containing a convenient number of tubes. Simply open one bag at a time, leaving the rest sealed, helping to prevent contamination.

PCR* plates, thin wall, for use with stripcaps **TUL-911-035R**

Catalogue No	Alt. No	Description	Volume, mL	Pack qty
TUL-910-010Q	3415	24 well plate	0.2	25

PCR* 8-tube strip, thin wall

Catalogue No	Alt. No	Description	Volume, mL	Pack qty
TUL-911-010K	3418	8 PCR strip tubes with flat caps, natural	0.2	120
TUL-911-035R	3418C	8 PCR strip tubes with dome caps, natural	0.2	120
TUL-911-040B	3418F	8 PCR strip tubes with flat caps, natural	0.2	120

PCR* tubes with caps, individual, pack of 1,000

Catalogue No	Alt. No	Description	Volume, mL
TUL-912-010E	3412	Flat caps, thin wall	0.2
TUL-912-011C	3412A	Flat cap, assorted colours	0.2
TUL-912-013V	3412G	Flat cap, green	0.2
TUL-912-015R	3412R	Flat cap, red	0.2
TUL-912-017N	3412B	Flat cap, blue	0.2
TUL-912-019J	3412P	Flat cap, purple	0.2
TUL-912-021W	3412Y	Flat cap, yellow	0.2
TUL-912-023S	3412O	Flat cap, orange	0.2
TUL-912-030V	3414	Dome caps, thin-wall	0.2
TUL-912-031T	3414A	Dome cap, assorted colours	0.2
TUL-912-033P	3414G	Dome cap, green	0.2
TUL-912-035L	3414R	Dome cap, red	0.2
TUL-912-037H	3414B	Dome cap, blue	0.2
TUL-912-039D	3414P	Dome cap, purple	0.2
TUL-912-041Q	3414Y	Dome cap, yellow	0.2
TUL-912-043M	3414O	Dome cap, orange	0.2
TUL-912-053J	3430G	Flat cap, green	0.5
TUL-912-055F	3430R	Flat cap, red	0.5
TUL-912-057B	3430B	Flat cap, blue	0.5
TUL-912-059U	3430P	Flat cap, purple	0.5
TUL-912-061K	3430Y	Flat cap, yellow	0.5
TUL-912-063G	3430O	Flat cap, orange	0.5
TUL-912-080G	3432	Dome cap, natural	0.5
TUL-912-081E	3432B	Dome caps, blue	0.5
TUL-912-083A	3432G	Dome cap, green	0.5
TUL-912-085T	3432O	Dome cap, orange	0.5
TUL-912-087P	3432P	Dome cap, purple	0.5
TUL-912-089L	3432R	Dome cap, red	0.5
TUL-912-091B	3432Y	Dome cap, yellow	0.5
TUL-912-050P	3430	Thin wall	0.5