



### Cassette labeller, Printmate™, Thermo Scientific Shandon

**Thermo**  
SCIENTIFIC

NEW

50

- Printmate prints easy-to-read, permanent black type (300dpi print is perfect for printing linear and 2-D barcodes as well as graphics and logos)
- The small size of the PrintMate means it can easily fit next to the grossing area and print fast enough to keep up when printing single cases; your printer allows you to decide the workflow that best fits your laboratory
- Prints cassettes with no reduction time for printing 2-D barcodes. Combines ruggedness with capacity to work reliably as a batch printer in high volume laboratories
- Dual track collection system collects up to 100 cassettes in organiser trays. The trays stacked and simplify matching cassettes to specimen containers at the grossing station
- LabWriter™ software prints large batches or single cassettes using standard or special numbering configurations. The software can use barcodes from requisitions to automatically print cassettes eliminating errors from the process
- The PrintMate is compatible with all LIS interfaces to any Thermo Scientific printer
- Connections includes: Serial, network, USB

Cassette marking system that permanently marks tissue cassettes with easy to read black type. Each unit comes with USB cable; power cable; 2, 6 or 12 hoppers (depending on model); one roll of tape (installed); basic software CD; and operator's guide. Complete cassette printing systems additionally include PC, monitor and digital imager.

#### Technical Specification - General

Dimensions [w x d x h], mm      480 x 400 x 230 (for Printmate 900, without cassette tubes)  
 480 x 400 x 900 (for Printmate 900, with cassette tubes)

Mass, kg      12

Please contact Fisher Scientific to discuss your EU compatible product requirements.

Catalogue No	Alt. No	Description
<b>SDA79500001</b>	A79500001	Printmate 150 cassette printer (2 hopper unit)
<b>SDA79500004</b>	A79500004	Printmate 150 cassette printing system
<b>SDA79500002</b>	A79500002	Printmate 450 cassette printer (6 hopper unit)
<b>SDA79500005</b>	A79500005	Printmate 450 cassette printing system
<b>SDA79500003</b>	A79500003	Printmate 900 cassette printer (12 hopper unit)
<b>SDA79500006</b>	A79500006	Printmate 900 cassette printing system

#### Accessories

Catalogue No	Alt. No	Description	Pack qty
<b>SDA79510013</b>	A79510013	Hot foil tape	5
<b>SDA79510014</b>	A79510014	Output collection system	1

### Cassette labeller, carousel-type, Thermo Scientific Shandon Carousel Cassette MicroWriter™

**Thermo**  
SCIENTIFIC

50

- Uses foil tape and a heated stylus to imprint characters that maintain legibility even through processing solutions
- Made to withstand the workload of busy laboratories, marking approximately 20,000 characters per roll of tape
- Printing can be inverted for reading in either direction
- Change cassette colour and style easily with fully automated dispensing and labelling process
- Easy to change hoppers provide the flexibility of allowing an unlimited number of colours and cassettes
- Imprints, letters, numbers, Roman numerals and special characters on the cassettes in four font sizes
- Ability to print barcodes on cassettes
- Prints on all 45° flat-faced cassettes
- Choice of collection systems makes organising cassettes a seamless operation in any laboratory environment

Included with unit: CD-ROM with unlicensed demo copy of Thermo Scientific Shandon MicroWriter™ software; Dispenser hoppers – six included; Pre-filled cassette tube sample pack; MicroWriter™ specific RS232 cable; Power cord; Operator manual.

